

MICROSCOPES

Compound Microscopes

Upright
Jenco



BC-211

BC-211 \$1,215.00

- Widefield 10X/20 eyepieces
- Binocular, 30° inclined, Siedentopf design viewing head with single diopter; 55-75 mm interpupillary adjustment
- DIN, plan achromat 4X, 10X, 40XR, and 100XR objectives on 5-place, inward facing nosepiece
- Ball-bearing mechanical stage is 20.9 cm(W) x 14 cm(D)
- Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
- 1.25 N.A. Abbe condenser with rack and pinion focusing, and iris diaphragm
- 20 watt halogen Koehler-type illuminator
- 27.8 cm(D) x 20 cm(W) x 40.9 cm(H)
- 5-year mechanical & optical/2-year electrical warranty



Meiji

MT5300H \$2,835.00

- Widefield 10X/22 eyepieces
- Trinocular, 30° inclined, Siedentopf design viewing head has single diopter; 53-75 mm interpupillary adjustment
- Infinity-corrected, U. Plan achromat 4X, 10X, 40XR, and 100XR objectives on 5-place, inward facing nosepiece
- Ball-bearing, ceramic coated mechanical stage is 19.1 cm(W) x 12.6 cm(D)
- Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
- 1.25 N.A. Abbe condenser has rack and pinion focusing, iris diaphragm, and filter holder
- 30 watt halogen Koehler-type illuminator
- 39 cm(D) x 21 cm(W) x 41 cm(H)
- Lifetime warranty

162-ASC



National

162-ASC \$747.00

- Widefield 10X/18 eyepieces; pointer standard in one eyepiece
- Binocular, 30° inclined, Jentsch design viewing head with dual diopters; 54-76 mm interpupillary adjustment
- DIN, super high contrast achromat 4X, 10X, 40XR, and 100XR objectives on 4-place, inward facing nosepiece
- Ball-bearing mechanical stage is 14 cm(W) x 13.5 cm(D); up-stop adjustment is provided
- Coaxial coarse and fine focus controls; coarse has tension adjustment
- 1.25 N.A. Abbe condenser with rack and pinion focusing, iris diaphragm, and swing-out filter holder
- 20 watt halogen illuminator
- 23.5 cm(D) x 18.5 cm(W) x 39.5 cm(H)
- Lifetime warranty

163-ASC \$956.00

165-ASC \$849.00



165-ASC

Inverted

Jenco



CP-2A1

CP-2A1 \$2,970.00

- Widefield 10X/20 eyepieces
- Trinocular, 30° inclined, Siedentopf design viewing head with single diopter; 55-75 mm interpupillary adjustment
- 6 DIN, LWD plan achromat objectives: brightfield 10X, 25X, and 40X; phase 10X, 25X, and 40X
- 5-place, side facing nosepiece
- Attachable mechanical stage includes four inserts
- Coaxial coarse and fine focus controls; coarse has tension and up-stop adjustments
- 1.25 N.A. Abbe condenser with rack and pinion focusing, iris diaphragm, and filter holder; working distance of 55 mm for phase contrast, or 125 mm if swung out for brightfield
- Centerable, dedicated phase annuli mounted on sliders
- 30 watt halogen Koehler-type illuminator
- 54.6 cm(D) x 22.9 cm(W) x 57.2 cm(H)
- 5-year mechanical & optical/2-year electrical warranty

BC-2A1 \$5,215.00

BC-364 \$9,716.00

Meiji Access www.meijitechno.com
and call us to place an order

Continued