


LAB, SURGICAL & DENTAL INSTRUMENTS



Cell Lifter/Scraper, polypropylene, Bag/10. 377560010\$15.50



Thumb Grip



Scissor Style

Ear Punches, 2 mm hole

Thumb Grip, 5 cm long, plated steel. 1538\$8.25
 Scissor Style, 11 cm long, stainless steel. 500077\$28.00

Forceps (Tweezers) & Tongs



Fine Point, Angled



Medium Point, Straight



Medium Point, Angled



Fine Point, Straight



Broad Point, Straight

Forceps, general

Fine Point, stainless steel, serrated handles

Straight, Serrated Tips

11.5 cm, guide pin, economy. 300-130\$2.60
 14 cm, Germany. 57.598\$26.60
 16.5 cm, black matte, Dumont. TWZ-103.70\$21.10
 20.5 cm, guide pin, economy. 300-141\$7.50
 25.4 cm, guide pin, economy. 300-143\$8.50

Straight, Smooth Tips, 20.5 cm, guide pin. 300-111\$8.50

Angled, serrated tips, 15.5 cm, guide pin. 3600\$4.00

Curved, serrated tips, 11.5 cm, guide pin. 300-131\$2.70

Medium-Fine Point, stainless steel, serrated handles

Curved, smooth tips, 15.2 cm. 300-011\$6.90

Medium Point, stainless steel, serrated tips & handles

Straight

11.5 cm, Economy. 2200\$2.25

11.5 cm, Premium, Germany. 540-035\$10.50

15.5 cm (6"), Economy. 2230\$3.15

15.5 cm (6"), Premium, Germany. 540-050\$12.60

20.5 cm (8"), Economy. 2240\$4.05

20.5 cm (8"), Premium, Germany. 540-055\$21.00

25.5 cm (10"), Economy. 2250\$4.80

25.5 cm (10"), Premium, Germany. 540-060\$31.50

30.5 cm (12"), economy. 2260\$5.85

Curved, 11.5 cm, guide pin, economy. 3150\$2.40

Continued

Student, stainless steel. 300-055\$1.70
 straight, blunt serrated tips, 12.5 cm
 Broad Point, stainless steel. EA-11\$1.00
 straight, serrated tips & handles, 13 cm
 Plastic, straight, serrated tips & handles, 13 cm. 97-3901\$0.45
 Corrosion and heat resistant ABS plastic is non-conductive,
 and non-contaminating. Will not scratch optics.



Styles: 5, 3c, 6, 7, 532, SS,
 AA, 3, AF, 2, OOB, 2A

Forceps, Micro, smooth tips

Economy, stainless steel

#2A Style, straight, 12 cm. M2A\$4.40

#3 Style, straight, non-magnetic, 12 cm. TWZ-736.16\$3.80

#5 Style, straight, non-magnetic, 11.5 cm. TWZ-736.22\$3.80

#7 Style, curved, non-magnetic, 11.5 cm. TWZ-736.28\$3.80

Economy Set, 12-pc. TWZ-700.12\$30.00

styles: 5, 3c, 6, 7, 532, SS, AA, 3, AF, 2, OOB, 2A



Style #3 (very sharp)



Style #5 (needle sharp)



Style #7 (sharp)

Dumont Hi-Tech, Inox, matte finish, Switzerland

#3 Style, straight, 12 cm. TWZ-301.16\$19.00

#5 Style, straight, 11 cm. TWZ-301.22\$24.10

#7 Style, curved, 11.5 cm. TWZ-301.28\$23.45

SS Style, straight, 14 cm. TWZ-301.56\$18.45

Dumont Hi-Tech, Titanium, non-magnetic, Switzerland

#3 Style, straight, 12 cm. TWZ-120.16\$28.00

Dumont Grippits, stainless steel, non-magnetic, Switzerland

#7 Style, serrated handles. TWZ-322.28\$19.56

Dumont Original, Inox, Switzerland

#3 Style, straight, 12 cm. TWZ-201.16\$26.16

#5 Style, straight, 11 cm. TWZ-201.22\$27.05

#7 Style, curved, 11.5 cm. TWZ-201.28\$29.50

Bishop Harmon



Cutting Forceps



Forceps, specialty, stainless steel

Adson Dressing, serrated tips and handles

Fine Tips, 12 cm. 17-555\$6.00

Medium Tips, 12 cm. 2360\$4.50

Adson Tissue, 1 x 2 teeth, 12 cm. 17-560\$6.00

Bishop Harmon, fine serrated tips, 9 cm. 4211\$18.00

Coverslipsee COVERSLLIP ACCESSORIES

Cutting, style 15A, Switzerland

Carbon Steel, 11 cm, Peer-Vigor. 57.512\$72.00

Stainless Steel, 11.5 cm, Dumont. TWZ-301.60\$39.85

Continued