




**Bone Cutters, stainless steel**

Liston

- 14 cm (5½") Long, 20 mm blades. **320-086**.....\$21.00
- 18.5 cm (7¼") Long, 25 mm blades. **1095**.....\$28.75
- Scissor Style, 10 cm (4") long, 18 mm blades. **320-025**.....\$8.75  
straight, strong tips



**Bone Nipper**, 10 cm (4") long, 6 mm cutting edge. **B286**.....\$12.00  
angled offset jaws


**Scalpels**




**Scalpels**

Disposable, plastic handle, sterile, Bx/10

- #10 Carbon Steel Blade. **9410**.....\$9.60
- #11 Carbon Steel Blade. **9411**.....\$9.60
- #15 Carbon Steel Blade. **9415**.....\$9.60
- #20 Carbon Steel Blade. **9420**.....\$9.60
- #22 Carbon Steel Blade. **9422**.....\$9.60



**Gel Scalpel**, Pk/10. **369700010** .....\$17.50  
clear polystyrene blade & handle



10 11 12 15 20 21 22 24 23

**Scalpel Blades**

Sterile

- Carbon Steel, Bx/100.....\$19.20
- #10 blade. **9310**                      #21 blade. **9321**
- #11 blade. **9311**                      #22 blade. **9322**
- #12 blade. **9312**                      #23 blade. **9323**
- #15 blade. **9315**                      #24 blade. **9324**
- #20 blade. **9320**

- Stainless Steel, Bx/100.....\$23.00
- #10 blade. **36-010S**                      #21 blade. **36-021S**
- #11 blade. **36-011S**                      #22 blade. **36-022S**
- #12 blade. **36-012S**                      #23 blade. **36-023S**
- #15 blade. **36-015S**                      #24 blade. **36-024S**
- #20 blade. **36-020S**

Stainless Steel, #60 blade, Bx/100. **36-060S** .....\$47.60  
Nonsterile, Swann-Morton, England

- Carbon Steel, Bx/100 .....\$22.00
- #10 blade. **KNF-245.10**
- #11 blade. **KNF-245.11**
- #12 blade. **KNF-245.12**
- #15 blade. **KNF-245.15**



**Scalpel Blade Remover**, holds 150-300 blades. **KNF-100.00**.....\$3.70




Stainless Steel, #3 Style



Stainless Steel, #4 Style



PVC Coated Handle, #3 Style



Plastic Handle, #5/6/8 Styles

**Scalpel Handles, reusable**

Stainless Steel, solid

- #3 Style, accepts blades #10-17
- Economy. **36-405** .....\$2.75
- Premium. **565-405** .....\$9.80
- #4 Style, accepts blades #18-36
- Economy. **36-410** .....\$2.75
- Premium. **565-410** .....\$11.20
- #7 Style, accepts blades #10-17
- Economy. **36-415** .....\$4.50
- Premium. **565-415** .....\$14.00

Stainless Steel, hollow

#8 Style, accepts blades #18-36. **1608** .....\$7.50

PVC Coated Stainless Steel

#3 Style, accepts blades #10-17. **KNF-225.00**.....\$5.35

Plastic, contoured to fit palm of hand

#5 Style, accepts blades #10-17. **KNF-221.00**.....\$5.00

#6 Style, accepts blades #18-36. **1605** .....\$5.00

#8 Style, accepts #60 blade, Bard-Parker™. **371080** .....\$32.40

Continued