



Scissors


George Tiemann & Co, GermanyAccess www.georgetiemann.com
and call us to place an order



Scissors, Micro, spring action, ambidextrous, stainless steel
Barraquer (butterfly), 6 cm long. **3737**\$34.50
10 mm angled-on-flat blades




Castroviejo, straight, sharp tips, economy
7 cm Long, 5 mm blades. **36205**\$36.00
Effective cutting edge is 3 mm; tip-to-tip is 3 mm.
8 cm Long, 7 mm blades. **5-254**\$36.00
Effective cutting edge is 6 mm; tip-to-tip is 4 mm.
8.5 cm Long, 5 mm blades. **5-855**\$36.00
Effective cutting edge is 4 mm; tip-to-tip is 3 mm.
11.5 cm Long, 17 mm blades. **M3720**\$26.00
Effective cutting edge is 12 mm.



Lifetime guarantee
against rust and defects.


Uniband™, sharp tips, Laschal Surgical, USA
Straight, 12 mm Blades. **LA-4**\$26.55
Straight, 20 mm Blades. **LA-1**\$26.55
Straight, 20 mm Serrated Blades. **LA-2**\$26.55
Curved, 12 mm Blades. **LA-4C**\$26.55
Curved, 20 mm Blades. **LA-1C**\$26.55
Featherlite™, Castroviejo, sharp tips, Laschal Surgical, USA
Greater tactile sensitivity than Uniband™ scissors.
Straight, 12 mm Blades. **MPF-4**\$80.80
Straight, 20 mm Blades. **MPF-1**\$80.80
Curved, 12 mm blades. **MPF-4C**\$80.80



Scissors, Dissecting, stainless steel
Straight, sharp tips, 11.5 cm (4½")
25 mm Blades, Economy. **12-410**\$5.05
25 mm Blades, Premium. **500-410**\$25.20
30 mm Blades, Economy. **310-040**\$4.50
35 mm Blades, Economy. **310-045**\$4.35
35 mm Blades, Premium. **510-045**\$21.00

Continued

Straight, Enamel-Coated, sharp tips, 11.5 cm (4½")
Economy. **MT30** *\$6.00
(*Patterns vary on these autoclavable tools.
Specify Leopard, Red Checker, or Zebra pattern.)




Curved, sharp tips, 11.5 cm (4½")
25 mm Blades, Economy. **310-035**\$4.35
25 mm Blades, Premium. **500-425**\$25.20
Curved, Enamel-Coated, sharp tips, 11.5 cm (4½")
Economy. **MT35** *\$6.00



(*Patterns vary on these autoclavable tools.
Specify Leopard, Red Checker, or Zebra pattern.)




Angled, 25 mm angled-on-side blades, economy. **310-066**\$6.00




Scissors, Operating, stainless steel
Straight, S/S (sharp/sharp) tips
14 cm (5½"), Economy. **12-055**\$5.15
14 cm (5½"), Premium. **510-055**\$23.80
15.5 cm (6"), Economy. **12-060**\$5.15
15.5 cm (6"), Premium. **510-060**\$23.80
16.5 cm (6½"), Economy. **12-065**\$6.75
16.5 cm (6½"), Premium. **510-065**\$28.00
Straight, S/B (sharp/blunt) tips
11.5 cm (4½"), Economy. **12-015**\$4.90
14 cm (5½"), Economy. **12-025**\$5.15
14 cm (5½"), Premium. **510-025**\$23.80
16.5 cm (6½"), Economy. **12-035**\$6.75
16.5 cm (6½"), Premium. **510-035**\$28.00
Straight, B/B (blunt/blunt) tips
11.5 cm (4½"), Economy. **12-075**\$4.90
14 cm (5½"), Economy. **12-085**\$5.15
14 cm (5½"), Premium. **510-085**\$23.80
16.5 cm (6½"), Economy. **12-095**\$6.75
16.5 cm (6½"), Premium. **510-095**\$28.00
Curved, S/S (sharp/sharp) tips
14 cm (5½"), Economy. **12-155**\$5.60
14 cm (5½"), Premium. **510-155**\$25.20
16.5 cm (6½"), Economy. **12-165**\$7.10
16.5 cm (6½"), Premium. **510-165**\$29.40
Curved, S/B (sharp/blunt) tips
11.5 cm (4½"), Economy. **12-115**\$4.90
14 cm (5½"), Economy. **12-125**\$5.60
14 cm (5½"), Premium. **510-125**\$25.20
16.5 cm (6½"), Economy. **12-135**\$7.10
16.5 cm (6½"), Premium. **510-135**\$29.40
Curved, B/B (blunt/blunt) tips
11.5 cm (4½"), Economy. **12-175**\$4.90
14 cm (5½"), Economy. **12-185**\$5.60
16.5 cm (6½"), Economy. **12-195**\$7.10
16.5 cm (6½"), Premium. **510-195**\$29.40

Continued


- Scissors, Metzenbaum, 15 cm (5¾"), stainless steel
 Straight, economy. 12-225 \$9.00
 Curved, economy. 12-230 \$9.35



- Scissors, Strabismus, 10.5 cm (4"), stainless steel
 Straight, economy. 26-105 \$6.25
 Curved, economy. 26-110 \$6.25



- Scissors, Student, 11.5 cm (4½"), steel wire, USA. 22 \$0.90
 straight, sharp/blunt tips



- Scissors, Bandage, stainless steel, economy
 Knowles, straight, sharp probe tips, 14 cm (5½"). 1070 \$6.25
 Lister, angled, blunt probe tips, 14 cm (5½"). 1050 \$5.00
 Can be used as enterotomy scissors.
 Lister, Enamel-Coated, angled, blunt probe tips, 14 cm (5½")
 14 cm (5½"). MT121 * \$6.75
 19 cm (7½"). MT120 * \$7.90
 (*Patterns vary on these autoclavable tools.
 Specify Leopard, Red Checker, or Zebra pattern.)



Cut through everything from wire to bone with these durable scissors. Bottom blade is serrated and has a footplate for safety. Blades are stainless steel & handles are plastic. Autoclavable to 120° C.

- Utility Shears, 18 cm (7¼"), economy. 310-341 \$5.25

Scrapers

Remove paraffin, ice, etc., safely from your laboratory equipment with these plastic scrapers.



- Scrapers, plastic
 5 cm (2") Wide Blade. 311-245 \$3.70
 10 cm (4") Wide Blade. 311-250 \$4.70

Spatulas & Scoops



- Spatulas, stainless steel blade, hardwood handle
 10 cm (4") Long Blade. 14-365B \$6.60
 15 cm (6") Long Blade. 14-365D \$8.60
 20 cm (8") Long Blade. 14-365E \$8.90



- Spatulas, Micro
 Square End/Spoon End, stainless steel
 19 cm (7½") Long. 320-411 \$4.15
 23 cm (9") Long. 320-414 \$4.25
 Round End/Tapered End. 320-370 \$4.65
 Round end: 8 x 51 mm; Tapered end: 6 x 51 mm
 nickel-stainless steel, 18 cm (7") long
 V-Shaped End/Round End. HS15906 \$4.95
 nickel-stainless steel, 15 cm (6") long



- Scoops, v-shaped
 Polypropylene, 22.5 cm (9"), Bel-Art, Pk/12. 378400000 \$21.55
 Stainless Steel, "reagent digger"
 23 cm (9") Long, Bel-Art, each. 367130000 \$20.70
 16.5 cm (6½") Long, each. 7-1835-1 \$2.25
 Handle, wooden, 13 cm long, each. 7-1836 \$2.85 ↓

Teasing Needles, Probes & Pin Vises



- Teasing Needles
 Wooden Handle, 14.5 cm (5¾")
 Angular. 97-1302 ea \$0.38 pk/12 3.80
 Straight. 97-1301 ea \$0.38 pk/12 3.80



- Stainless Steel, straight, 14 cm. M28 \$2.20
 Probe/Seeker, stainless steel, 15.5 cm. 320-216A \$4.50
 One end tapered to a point; one end curved with a rounded tip.

Continued