



Probe with Eye, flexible, 6". 5850\$2.95




Hold needles, probes, wires, drill bits, etc.

Pin Vise, Swivel-Head, 0-3.3 mm opening. PIN-220.00\$7.95

Two double-ended collets are provided; one in screw chuck, one inside handle.

Pin Vise, Aluminum, 0-3.0 mm opening. PIN-222.00.....\$4.10

Trocar



Trocar, 4-in-1 set, 13 cm, stainless steel. 1596\$17.50

Accessories/Maintenance/Storage



Cutting Boards, white polyethylene


10 x 6 x 1/2" (25.4 x 15.2 x 1.3 cm). JRC4319\$6.95

18 x 12 x 1/2" (45.7 x 30.5 x 1.3 cm). AAACBP12.....\$13.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1 1/2 x 3" opening. DEM-345.00\$62.95



Dissecting Cases

6-Loop with snap fastener, black leatherette. 7-370.....\$1.10

12-Loop with snap fastener, black leatherette. 7-371.....\$1.50

13-Loop with zipper, brown leatherette. 36900.....\$7.80

Features red velveteen interior, protector flap, and ID pouch.

Dissecting Kit, Basic. 38501.....\$5.30


Contents: Dissecting Case, 6-loop, black vinyl (7-370); Forceps, broad point, stainless steel, 13 cm (EA-11); Pipet, plastic transfer, 5 mL capacity (PTP-02); Pipet, plastic transfer, 7 mL capacity (PTP-01); Scalpel, disposable, #10 carbon steel blade (9410); Scissors, steel wire, 11.5 cm (22); Teasing needle, wooden handle, straight (97-1301); Teasing needle, wooden handle, angular (97-1302); and Ruler, white vinyl, 15 cm/6".

Dissecting Pans, aluminum with wax

Standard, 11 x 7 x 1 1/2". 7-350.....\$13.00

Large, 13 1/8 x 9 3/8 x 2 1/4". 7-355.....\$19.50

Replacement Wax, black, 5 lb Pk. 7-354-1\$16.20



Instrument Sterilizing Trays, plastic, autoclavable, PST, USA

Base, lid & mat (B/L/M) included.


6 x 2 1/2 x 3/4" (15.2 x 6.4 x 1.9 cm). 2600A\$20.00

(500252)

6 1/2 x 4 x 3/4" (16.2 x 10.2 x 1.9 cm). 4600A\$29.00

10 x 6 x 3/4" (25.4 x 15.2 x 1.9 cm). 6100A\$60.00

(500254)



Base, insert tray, lid & 2 mats (B/I/L/2M) included.

15 x 10 x 1 1/2" (38.1 x 25.4 x 3.8 cm). 1520A\$398.00



Instrument Storage Box, 10 x 4 x 2", stainless steel. 6830\$19.50

Tools




Cutter, diagonal, 5". PLR-766.00 \$10.95




Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1½ x 3" opening. DEM-345.00 \$62.95

Gauges



Vernier Caliper



Dividers

Vernier Caliper, digital, 3-button. GAU-178.00 \$39.95
 Powered by a 1.5 V cell and protected in a foam-lined, plastic case. Measures in millimeters (0-100 mm with 0.01 mm precision) and inches (0-4").

Replacement Battery. 357 \$4.00

Dividers (compasses), replaceable tips, Germany
 3" Long. DIV-320.00 \$11.90
 6" Long. DIV-326.00 \$16.95
 Replacement Tips, Pair. DIV-320.05 \$1.85



Reproducibility enhanced by rapid stop. Provided with padded, plastic case.

Micrometer, 0-25 mm in 0.01 mm div, Germany. 601A \$39.00



Measure inside and outside diameters of glass and plastic tubing, connecting pieces, rods, test tubes, etc. between 1.5 and 16 mm.

Tubing Sizer, metric diameter gauge, PP. 133900000 \$25.80

Hex L-Wrenches, (Allen), metric

- 0.71 mm, short-arm, each. 78001 \$0.75
- 0.89 mm, short-arm, each. 78002 \$0.75
- 1.27 mm, short-arm, each. 78003 \$0.75
- 1.5 mm, short-arm, each. 78004 \$0.75
- 2 mm, short-arm or long-arm (*specify*), each. 78005 \$0.75
- 2.5 mm, short-arm, each. 78006 \$0.75
- 3 mm, short-arm or long-arm (*specify*), each. 78007 \$0.75

Hex L-Wrench Sets

- 7-Pc, short-arm. 78000 \$5.00
- 0.71, 0.89, 1.27, 1.5, 2, 2.5 and 3 mm. Includes case.




Ball end can be inserted up to 20° off-axis and not bind. Includes storage case.

- 13-Pc, ball end, long-arm, chrome vanadium
- Metric. 96416 \$15.90
- 1.5, 2, 2.5, 3, 3.5, 4, 4.5, 5, 5.5, 6, 8, and 10 mm
- SAE. 94680 \$15.90
- 0.050, 1/16, 5/64, 3/32, 7/64, 1/8, 9/64, 5/32, 3/16, 7/32, 1/4, 5/16, and 3/8"



Inspection Mirror, 9", stainless steel. M31 \$2.50



Needle File Set, 12-pc, 10 cm long. FIL-990.00 \$15.40
 Needle File Handle, anodized aluminum. HAN-739.00 \$11.50