


Fused Quartz, USA

0.15-0.25 mm Thick, NSG Type OZ®	
1" Square, 1 each. QCS1X1	\$19.00
1" Circle, 1 each. QCS1X0	\$19.00
0.16 ± 0.01 mm Thick, A1	
25 mm Square, 1 each. R525025	\$32.40
25 mm Circle, 1 each. R525000	\$23.40

Fused Silica, USA

0.15-0.25 mm Thick, Corning 7980	
1" Square, 1 each. FSCS1X1	\$21.50
1" Circle, 1 each. FSCS1X0	\$21.50
0.16 ± 0.01 mm Thick, S1-UV	
25 mm Square, 1 each. R425025	\$32.00
25 mm Circle, 1 each. R425000	\$28.80

**Accessories**



	<u>bk</u>	<u>48 bk</u>
Bibulous Paper, 10 x 15 cm, Bk/50 sheets. PB46L	\$3.00	130.00

Forceps



Minimize contamination during slide & coverslip handling.

Kirkbridge Slide



Hemostat, polypropylene, 12.5 cm. HEM-300.00	\$1.25
Will not scratch; corrosion and heat resistant.	
Kirkbridge Slide, nickel-plated wire, 12.5 cm long. 7-G36	\$1.95



Spade Tip



Square Point Tip

Spade Tips

Angled (Kuehne)	
Stainless Steel, 11 cm (4¼"). 3270	\$4.15
Teflon® Coated, 11 cm (4¼"), Bel-Art. 379420000	\$16.65
Straight, stainless steel. 57.926	\$4.90

Square Point Tips, straight

Stainless Steel, 12 cm (4¾"). 1244	\$2.25
Teflon® Coated, 11 cm (4¼"), Bel-Art. 379430000	\$18.70

**Immersion Oils, Mounting Media & Refractive Index Liquid**



Cargille immersion oils A, B, and NVH can be mixed with each other for custom viscosities. USA

**Type A, low viscosity (150 cSt), 16482. Short W.D. optics**

7.5 mL in LDPE dispensing bottle. 18101	\$3.40
20 mL in HDPE dispensing bottle. 18102	\$8.30
30 mL in amber bottle w/applicator rod. 18109	\$6.60
120 mL in amber bottle. 18103	\$15.30
480 mL in amber bottle. 18104	\$33.20

**Type B, high viscosity (1250 cSt), 16484. Long W.D. optics**

7.5 mL in LDPE dispensing bottle. 18105	\$3.40
20 mL in HDPE dispensing bottle. 18106	\$8.30
30 mL in amber bottle w/applicator rod. 18110	\$6.60
120 mL in amber bottle. 18107	\$15.30
480 mL in amber bottle. 18108	\$33.20

**Type NVH, very high viscosity (21,000 cSt), 16485**

120 mL in amber bottle. 18112	\$20.35
480 mL in amber bottle. 18113	\$48.50

**Type LDF (formerly DF), low fluorescence (500 cSt), 16241**

7.5 mL in LDPE dispensing bottle. 18115	\$5.35
20 mL in HDPE dispensing bottle. 18116	\$9.40
120 mL in amber bottle. 18117	\$25.45
480 mL in amber bottle. 18118	\$56.50

**Type FF, fluorescence free (170 cSt), 16212**

7.5 mL in LDPE dispensing bottle. 18121	\$4.85
20 mL in HDPE dispensing bottle. 18122	\$9.25
120 mL in amber bottle. 18123	\$23.05
480 mL in amber bottle. 18124	\$54.10

**Type 37, 37°C (1250 cSt @ 37°C), 16237**


20 mL in Dispensing Bottle. 18125	\$9.00
120 mL in amber bottle. 18126	\$20.35
480 mL in amber bottle. 18127	\$48.50

**Non-Volatile Water Substitute (refractive index liquid)**

Substitute this liquid for water to minimize evaporation when using water immersion objectives. Refractive index is 1.330 at 23°C; viscosity 15 cSt.

7.5 mL in amber bottle w/applicator rod. 1803X	\$45.50
30 mL in amber bottle w/applicator rod. 1803Y	\$115.00

A long nozzle, textured sides, and slim profile make this chemically inert, high-density polyethylene (HDPE) container ideal for dispensing immersion oil.



Dispensing Bottle, HDPE, 20 mL. 18145	\$6.10
---------------------------------------	--------



Store immersion oil or mounting media in this borosilicate glass bottle with an applicator rod and ground glass cap.

Balsam Bottle, glass, 45 mL (empty), USA. 616	\$16.00
---	---------

Continued

# GENERAL LAB SUPPLIES & EQUIPMENT

## Applicators, 6" long



Cotton Tipped, wooden shaft, nonsterile.

- Bg/100. 95-8702A ..... \$1.25
- Bx/1000. 95-8702B..... \$10.20
- Plain, wood, Bg/1000. 95-8701 ..... \$6.00



Bag Sealer, electric, 12" (30.5 cm) Wide. HS28272..... \$180.00 ↓

SCL-292.50




## Balances, portable

- ★ 1000 g Capacity x 1 g. XH-400 ..... \$31.50  
Operates on 2 AA batteries; 5 3/4" white plastic platform.
- 1000 g Capacity x 0.1 g. SCL-292.50 ..... \$80.00  
Operates on 4 AA batteries or AC adapter; 5 3/4" SS platform.

## Beakers

Glass (borosilicate), Griffin low form, single scale, economy

- 50 mL, each. 7-440050..... \$2.15
- 100 mL, each. 7-440100..... \$2.35
- 250 mL, each. 7-440250..... \$3.00
- 600 mL, each. 7-440600..... \$4.50
- 1000 mL, each. 7-441000..... \$8.40
- 2000 mL, each. 7-442000..... \$17.25
- 4000 mL, each. 7-444000..... \$32.90



Economy



Premium



Stainless Steel

Polypropylene

- | Economy, three pouring spouts, Tri-Stir® | ea     | cs/100  |
|--|--------|---------|
| 50 mL, 5 mL graduations. 80092.....      | —      | \$25.75 |
| 100 mL, 10 mL graduations. 80093.....    | —      | \$30.20 |
| 250 mL, 20 mL graduations. 80094.....    | \$0.85 | 40.30   |
| 1000 mL, 50 mL graduations. 80098.....   | \$1.75 | 79.60   |
- Premium, Griffin low form, single scale, Bel-Art
- 100 mL, 25 mL graduations, each. 262110000 ..... \$2.60
  - 250 mL, 50 mL graduations, each. 262130000 ..... \$4.60
  - 600 mL, 100 mL graduations, each. 262150000 ..... \$7.30
  - 1000 mL, 200 mL graduations, each. 262160000 ..... \$7.40
  - 2000 mL, 500 mL graduations, each. 262170000 ..... \$17.10
  - 4000 mL, 1000 mL graduations, each. 262190000 ..... \$33.00

Stainless Steel, griffin style, nongraduated, full capacity mark

- 125 mL, USA, each. 125B ..... \$34.90
- 250 mL, USA, each. 250B ..... \$35.25
- 600 mL, USA, each. 600B ..... \$40.50
- 1200 mL, USA, each. 1200B ..... \$50.55



Bench Liner, white, 50' x 20" roll, Labmat™. 246750000..... \$51.45



These LDPE bottles deliver controlled drops from standard 0.030" tip opening. Sealer cap snaps on easily after use.




## Bottles

Dispensing, squeeze bottles

- LDPE (low-density polyethylene), Yorker cap, Bel-Art
- 30 mL (1 oz), Bg/12. 116370001 ..... \$14.90
- 60 mL (2 oz), Bg/12. 116370002..... \$17.35
- 125 mL (4 oz), Bg/12. 116370004..... \$28.95
- 185 mL (6 oz), Bg/12. 116370006..... \$28.95

A long nozzle, textured sides, and slim profile make this chemically inert, HDPE container ideal for dispensing immersion oil.



- HDPE (high-density polyethylene)
- 20 mL, each. 18145..... \$5.80

★ Amber glass protects light-sensitive liquids. Screw-cap dropper features a black rubber bulb and glass pipet.



## Dropping

- Amber Glass, Boston Round
- 30 mL (1 oz), Bx/12. 7-405-3..... \$17.10
- 60 mL (2 oz), Bx/12. 7-405-4..... \$27.20



Glass dropper with vinyl bulb fits into polypropylene stopper to form a tight seal.



- Clear Glass, round, Polystop™, Bel-Art
- 50 mL, each. 116040060..... \$18.45
- 100 mL each. 116040125..... \$20.30

Continued

Continued