



Slide Sleeves, 10 gauge vinyl
 4 x 5½" Card Holder Plus 4 Slide Pockets. **SL-3-264**\$0.45
 File Box, 4 x 6 x 11½", fiberboard. **40589**\$16.80
 5 x 8" Card Holder Plus 5 Slide Pockets. **SL-3-108**\$0.55
 File Box, 6 x 9 x 11½", fiberboard. **40591**\$26.25
 8½ x 11" Page Holds 18 Slides. **SL-3-182**\$0.75



Slide Folders

NEW 12-Slide, holds 75 x 50 mm (3 x 2") slides ea cs/72
 Cardboard, black tape, USA. **SF12-BLACK**\$6.50 406.00
 20-Slide, holds 75 x 25 mm (3 x 1") slides ea cs/72
 Cardboard, USA\$6.50 406.00
 Black Tape. **SF20-BLACK**\$6.50 406.00
 Canary Yellow Tape. **SF20-YELLOW**\$6.50 406.00
 Fire Engine Red Tape. **SF20-RED**\$6.50 406.00
 Forest Green Tape. **SF20-GREEN**\$6.50 406.00
 Navy Blue Tape. **SF20-BLUE**\$6.50 406.00
Folders made in the USA resist heat-induced warping
 Cardboard, Black Tape, China, each. **HS9902**\$4.65
 (12-587-10)
 Plastic, gray tray w/2 clear doors, each. **62680**\$9.15
 Transparent covers hold visible slides securely in place,
 but can fold flat against back of tray when opened.
 Slides are slightly elevated for easy retrieval.



Slide Holders, perfect for slide drying or staining
 12-Slide (horizontally positioned), polystyrene. **445570000**\$6.60
 52-Slide (approx.), 12 x 4", HDPE, USA. **L500-C**\$27.20
 91-Slide (approx.), 12 x 7", HDPE, USA. **L500-D**\$36.35



Slide Mailers

Cardboard, top corner notched	<u>ea</u>	<u>cs/2000</u>
1-Slide. MSL-01	\$0.13	214.00 ↓
2-Slide. MSL-02	\$0.21	176.00 ↓
Plastic	<u>ea</u>	<u>cs/1000</u>
1-Slide. 880052	\$0.13	102.00
2-Slide. 880041	\$0.30	116.00
4-Slide, screw-on cap. 950240	—	90.00
5-Slide, end opening, USA. 62706	\$0.48	80.00

62706 mailers can also be used as staining jars.
 (1258717B)

Slide Pens, Markers & Scribes

Pens



Pigment Liner, 0.5 mm stroke, Staedtler	<u>ea</u>	<u>bx/10</u>
Black. 308-05-9	\$3.50	\$29.00

Write a 0.4 mm stroke on glass slides with this permanent marker.
 Ideal for cytologists and pathologists.



Pilot®, permanent, ultra fine tip	<u>ea</u>	<u>bx/12</u>
Black. 44102	\$1.90	19.60 ↓
Blue. 44103	\$1.90	19.60 ↓
Green. 44105	\$1.90	19.60 ↓
Red. 44104	\$1.90	19.60 ↓




Confine immunohistochemical reagents within a light blue, hydrophobic barrier, or separate multiple specimens on a slide/dish. Resists alcohol and acetone, but removes with xylene.

Super^{HT} PAP (Peroxidase Anti-Peroxidase), hydrophobic
 Regular Tip (4 mm tapered nib). **K-600**\$59.00
 (71310, 195505)
 Fine Tip (2.5 mm tapered nib). **K-650**\$44.00
 (71312, 195506)

Continued

Markers

Ink resists xylene, alcohol, acetone and formalin.
Remove with acetic acid, chloroform, or ethanol/xylene mixture.



Securline® Marker II®, permanent, PDC, USA ea bx/10
Black, 0.6 mm tip. 1450-20-PDC \$2.50 20.30
(14-905-30, 72107)





Set includes 8 water-resistant colors: berry, black, lime, navy, orange, purple, sky blue, and turquoise.

Sharpie Dual-Tip Set, fine & ultra-fine tips, Pk/8. 33861 \$19.95




Sharpie Retractable, permanent

Fine Tip		<u>ea</u>	<u>bx/12</u>
Black. 36701		\$2.45	26.70
Blue. 36703		\$2.45	26.70
Red. 36702		\$2.45	26.70
Ultra-Fine Tip		<u>ea</u>	<u>bx/12</u>
Black. 1735790		\$2.45	26.70
Blue. 1735792		\$2.45	26.70
Red. 1735791		\$2.45	26.70



Pen Holder, Captor™, Bel-Art, Pk/3. 133730000 \$17.20



Scribes, retractable

Carbide-Tipped (resharpenable). RSC	\$15.30
Replacement Carbide Tip. RCT-1	\$6.20
Diamond-Tipped. RSD	\$19.10
Replacement Diamond Tip. RDT-1	\$12.90

Aluminum



Plastic

Slide Trays

Aluminum, Eberbach, USA. 2870	\$20.60
Plastic, white, 20-slide, Bel-Art. 445810000	\$10.35
Wood, 50-slide	\$42.00



Slide Warmers


Small, 26 x 18 cm (10 x 7"). XH-2002	\$225.00
Large, 63 x 20 cm (25 x 8"). XH-2001	\$299.00

Staining

Glass staining jar with ground glass lid holds 8 coverslips, 18-22 mm, back-to-back. Jar is 45 mm wide and tall.



Coverslip Staining Jar. 200-011 \$11.35



Staining/Moisture Chambers

10-Slide, polystyrene, USA	<u>ea</u>	<u>pk/3</u>
Clear. 240-9000-010	\$48.20	135.00
Black. 240-9020-Z10	\$50.10	141.00



20-Slide, Simport

Clear Lid (PETG), black base (ABS), each. M920-1	\$159.00 ↓
Black Lid (ABS), black base (ABS), each. M920-2	\$159.00 ↓



Staining Dishes

Glass (flint)	
20-Slide (back-to-back), holds slides horizontally, USA	
Set (grooved dish & lid). 672	\$12.10 ↓
Bottom only. 672-BTM	\$8.95
Lid only. 672-LID	\$4.95



Set includes 8 water-resistant colors: berry, black, lime, navy, orange, purple, sky blue, and turquoise.

Markers

Sharpie Dual-Tip Set, fine & ultra-fine tips, Pk/8. 33861 \$19.95 ↓




Sharpie Retractable, permanent

	ea	bx/12
Fine Tip		
Black. 36701	\$2.45	26.70
Blue. 36703	\$2.45	26.70
Red. 36702	\$2.45	26.70
Ultra-Fine Tip		
Black. 1735790	\$2.45	26.70
Blue. 1735792	\$2.45	26.70
Red. 1735791	\$2.45	26.70



Micromanipulation Systems



Capillaries, custom-tooled, Sutter, (6 pk min) ... Access www.sutter.com and call us to place an order



Microforge, Staelting Access www.staeltingco.com/physio and call us to place an order




Micromanipulators, Staelting Access www.staeltingco.com/physio and call us to place an order



Micropipette Pullers, Staelting Access www.staeltingco.com/physio and call us to place an order

Mortar and Pestle Sets, glass, USA

59 mL (2 oz) Capacity. 826-2	\$14.05
118 mL (4 oz) Capacity. 826-3	\$14.75
237 mL (8 oz) Capacity. 826-4	\$17.50



Dimensions: 17" x 11.5" x 12"(H).

Organizer, benchtop, TrippNT

Corner Glove Box & Wiper Holder. 50130 \$80.00



Petri Dishes, polystyrene, Cs/500, sterile sleeve/20, Simport

35 x 10 mm, vented lid, stacking ring. D210-15	\$123.00 ↓
60 X 15 mm, stacking ring. D210-13	\$97.00 ↓
100 X 20 mm, stacking ring. D210-7	\$164.00



Holds (5) 100 mm diameter, or (15) 60 mm diameter petri dishes. Dimensions: 25.1 x 23.7 x 3.5 cm.

Petri Dish Incubation Tray, polypropylene, stackable, Bel-Art

Pk/3. 189831000 \$46.85



Store (54) 100 mm petri dishes in this clear, 6-column rack. 34(L) x 20(W) x 23 cm(H).

Petri Dish Rack, 100 mm, acrylic, Bel-Art. 189850100 \$106.00

Pipetting



Dropper Bulbs, amber latex

1 mL Capacity, 43 mm long, Pk/12. 20621A	\$5.60
1 mL Capacity, 43 mm long, Pk/72. HS20622A	\$29.25
2 mL Capacity, 52 mm long, Pk/12. 20621B	\$5.60
2 mL Capacity, 52 mm long, Pk/72. HS20622B	\$29.25

Continued