



20-Slide, holds slides horizontally	
Set (dish, rack & lid). <b>SDS-20</b> .....	\$25.10
Dish & lid. <b>SD-20</b> .....	\$15.30
Rack only, stainless steel. <b>SR-20</b> .....	\$12.60
30-Slide, holds slides vertically	
Set (dish, rack & lid). <b>SDS-30</b> .....	\$34.30
Dish & lid. <b>SD-30</b> .....	\$17.10
Rack only, stainless steel. <b>SR-30</b> .....	\$21.00
60-Slide, holds slides vertically	
Set (dish, rack & lid). <b>SDS-60</b> .....	\$49.40
Dish & lid. <b>SD-60</b> .....	\$27.00
Rack only, stainless steel. <b>SR-60</b> .....	\$27.90



Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).

Plastic (polyoxymethylene [POM]), holds 25 slides  
Set (dish, rack & lid). **0990-99**.....\$49.40



674



Bel-Art  
442081000

675



**Staining Jars**, hold slides (up to 76 x 26 mm) vertically

Coplin, holds 10 slides back-to-back		<u>ea</u>	<u>cs/12</u>
Glass w/Ground Glass Cover, USA. <b>673</b> .....	\$10.70	↓	109.00
Replacement Glass Cover. <b>673-LID</b> .....	\$4.15	—	
		<u>ea</u>	<u>cs/12</u>
Glass w/Screw-On Phenolic Cap, USA. <b>674</b> .....	\$10.70	↓	109.00
Replacement Phenolic Cap, black. <b>674-LID</b>	\$4.25	—	
		<u>ea</u>	<u>cs/12</u>
Polypropylene w/Screw-On Cap. <b>442081000</b> .....	\$7.45		75.35
Plastic jar is unbreakable, microwave proof, stain resistant, and chemically inert.			
Wide Top		<u>ea</u>	<u>cs/12</u>
Glass, holds 16 slides back-to-back. <b>675</b> .....	\$14.40	↓	147.00
grooved jar w/square lid, USA			
Lid only. <b>675-LID</b> .....	\$4.50	—	




Plastic (POM)

Move slides, without forceps, between jars linked together in a workstation. Swing out the rack's integral lid for easy accessibility to slides. Stable, wide-based jar requires only 80 mL of reagent. Polyoxymethylene (POM) resists alcohol and xylene; can be microwaved or frozen.

12-Slide, EasyDip™, Simport	<u>ea</u>	<u>cs/6</u>
Jar only, Blue. <b>M900-12B</b> .....	\$10.00	55.00
Jar only, Green. <b>M900-12G</b> .....	\$10.00	55.00
Jar only, Pink. <b>M900-12P</b> .....	\$10.00	55.00
Jar only, White. <b>M900-12W</b> .....	\$10.00	55.00
Jar only, Yellow. <b>M900-12Y</b> .....	\$10.00	55.00
Rack/Lid only, gray. <b>M905-12DGY</b> .....	\$11.10	61.20
Staining Jar Holder, 6 jar capacity. <b>M906</b> .....	\$61.85	
stainless steel		



1030



1042

25-Slide, light green	
Jar only. <b>1030</b> .....	\$11.35
Rack only	
Metal Handle. <b>1040</b> .....	\$15.20
Plastic Handle (microwaveable). <b>1042</b> .....	\$14.70


**Tissue (Embedding) & Biopsy Cassettes**



<b>Base Molds</b> , disposable, Simport	
15 x 15 x 5 mm, Cs/1000. <b>M475-2</b> .....	\$122.00
24 x 24 x 5 mm, Cs/1000. <b>M475-3</b> .....	\$122.00
30 x 24 x 5 mm, Cs/1000. <b>M475-4</b> .....	\$122.00
37 x 24 x 5 mm, Cs/1000. <b>M475-5</b> .....	\$122.00

**Biopsy Bags** (Tea Bags), paper, 2¼ x 2¼", Pk/500. **56100** ..... \$34.00

Made of 60 PPI cellular polyester urethane foam.



<b>Biopsy Pads</b>	
Black Foam, 30.2 x 25.4 mm, Bg/1000. <b>339014</b> .....	\$38.85
Blue Foam, 30.2 x 25.4 x 2 mm, Cs/10,000. <b>M476-1</b> .....	\$174.00

Continued