


Molded from a high-density, solvent resistant polymer. Textured writing surface on front and sides.

Cassettes, Detached Lids (packed separately)

- Biopsy, small square pores, economy, Pk/500.....\$39.50
 - Blue. **BC-0109B** White. **BC-0109W**
 - Green. **BC-0109G** Yellow. **BC-0109Y**
- Tissue, slotted pores
 - Economy, Pk/500.....\$39.50
 - Blue. **EC-0109B** Pink. **EC-0109P**
 - Green. **EC-0109G** White. **EC-0109W**
 - Lilac. **EC-0109L** Yellow. **EC-0109Y**
 - Uni-Cassettes™, premium, USA, Pk/500.....\$57.50
 - Aqua. **6059** Lime. **6062**
 - Blue. **6054** Peach. **6064**
 - Burnt Orange. **6056** Pink. **6052**
 - Chartreuse. **6061** Red. **6060**
 - Gold. **6058** Rose. **6063**
 - Green. **6053** Tan. **6055**
 - Lavender. **6065** White. **6050**
 - Lilac. **6057** Yellow. **6051**



Histosette I




Micromesh

Cassettes, Attached Lids, acetyl polymer, premium

- Biopsy
 - Histosette™ I, 1 mm square pores, Pk/1500 ...\$131.00 1 pk 3 pk 354.00
 - Aqua. **M491-12** Peach. **M491-7**
 - Blue. **M491-6** Pink. **M491-3**
 - Green. **M491-4** Tan. **M491-8**
 - Lilac. **M491-10** White. **M491-2**
 - Orange. **M491-11** Yellow. **M491-5**
- Micromesh™, 0.38mm pores, Pk/10001 pk 3 pk \$238.00 ↓ 637.00↓
 - Aqua. **M507-12** Peach. **M507-7**
 - Blue. **M507-6** Pink. **M507-3**
 - Green. **M507-4** Tan. **M507-8**
 - Lilac. **M507-10** White. **M507-2**
 - Orange. **M507-11** Yellow. **M507-5**




Histosette I



Macrosette

- Tissue, slotted pores
 - Histosette™ I, snap-apart lids, Pk/1500..... 1 pk 3 pk \$108.00 290.00
 - Aqua. **M490-12** Peach. **M490-7**
 - Blue. **M490-6** Pink. **M490-3**
 - Green. **M490-4** Tan. **M490-8**
 - Lilac. **M490-10** White. **M490-2**
 - Orange. **M490-11** Yellow. **M490-5**
- Uniset™, hinged lids, Pk/1500 1 pk 3 pk \$116.00 312.00
 (Alternative to Sakura's Uni-Cassette®)
 - Aqua. **M505-12** Peach. **M505-7**
 - Blue. **M505-6** Pink. **M505-3**
 - Green. **M505-4** Tan. **M505-8**
 - Lilac. **M505-10** White. **M505-2**
 - Orange. **M505-11** Yellow. **M505-5**
- Macrosette™, 2X normal depth, white, Pk/750. **M512** \$173.00



M480-4

Cassettes, without Lids (compatible with all available metal lids)

- Tissue, acetyl polymer, premium
 - Capsettes™, USA, Pk/500.....\$39.60
 - Aqua. **6039** Lime. **6042**
 - Blue. **6034** Peach. **6044**
 - Burnt Orange. **6036** Pink. **6032**
 - Chartreuse. **6041** Red. **6040**
 - Gold. **6038** Rose. **6043**
 - Green. **6033** Tan. **6035**
 - Lavender. **6045** White. **6030**
 - Lilac. **6037** Yellow. **6031**
 - Simport, Pk/1500
 - Green, **M480-4**\$89.50



Store 1500 cassettes or 990 embedding rings in 6 stackable drawers. High-impact polystyrene resists temperatures from -80°C to +80°C. Drawer dimensions: 40.5 x 23 x 5.1 cm.

- Cassette Drawers, polystyrene, Pk/6. **M495-6** \$147.00

- Freeze Spray, 10 oz Can. **403A-285G** \$11.00



Continued

Continued