




Custom blend of highly purified paraffin and plastic polymers formulated for use as both an infiltration and embedding medium. Melting point: 58-63°C.

Paraffin, infiltration & embedding, 1 kg/bag 

Bag. SF007-1.....\$9.40

Cs/8 Bags. SF007-8 .....\$68.00

Paraffin Scrapers - see Scrapers in LAB, SURGICAL, & DENTAL INSTRUMENTS section



Paraffin Block Mailers, 4-block, polystyrene, Pk/10. 48315.....\$21.50

Tissue Block Mailing Sleeve, Block Tunnel™ .....\$2.25

**Equipment**



Paraffin Dispenser, 2.5 gallon. XH-4002.....\$725.00 ↓



Tissue Floatation Bath. XH-1001.....\$184.00