

**Desiccators**



- Jar, glass w/porcelain sample plate, economy  
 Glazed porcelain plate has 3.5-4.5 cm center hole.
- 21 cm (8.25") ID. **7-200-001** ..... \$44.75
  - 24 cm (10") ID. **7-200-002** ..... \$70.45
- Vacuum, polycarbonate top/polypropylene bottom, Bel-Art
- 140 mm Plate Diameter. **420100000** ..... \$65.60
  - 190 mm Plate Diameter. **420200000** ..... \$87.45
  - 230 mm Plate Diameter. **420250000** ..... \$119.00

**Bel-Art Cabinets**




- Cabinets, acrylic
- Clear, 30.5 x 30.5 x 30.5 cm. **420660000** ..... \$487.00
  - Bronze, 30.5 x 30.5 x 30.5 cm. **420660001** ..... \$559.00



- Dippers**, 500 mL polypropylene cup with HDPE handle, Bel-Art
- 3' Handle. **367840016** ..... \$42.75
  - 6' Handle. **367850016** ..... \$55.75
  - Replacement Cup, screw-on. **367860016** ..... \$9.95

- Dish**, specimen, glass, stackable, USA. **600** ..... \$9.15 ↓  
 120 mm(OD) x 50 mm



- Dry Baths**, digital, Benchmark
- 1-Block Capacity. **BSH1001** ..... \$339.00
  - 2-Block Capacity. **BSH1002** ..... \$427.00
  - 4-Block Capacity. **BSH1004** ..... \$759.00
- Blocks & bath beads sold separately – please specify configuration*
- Blocks, aluminum
- Holds (24) 1.5 mL Tubes. **BSW1500** ..... \$93.00
  - Holds (24) 1.5/2.0 mL Tubes. **BSW1520** ..... \$93.00
  - Holds (24) 0.5 mL Tubes. **BSW05** ..... \$93.00
  - Holds (12) 15 mL Tubes. **BSW15** ..... \$93.00
  - Holds (5) 50 mL Tubes. **BSW50** ..... \$93.00
- Lab Armor™ Beads, metallic thermal media, USA
- 0.75 L (replaces two blocks in a dry bath). **42370-750** ..... \$98.00
  - 2 L. **42370-002** ..... \$191.00

**Drying Racks**



- Non-Electric
- Bel-Art. **188050000** ..... \$89.00
- Dry labware on this epoxy-coated steel wire rack with polystyrene tray. Holds up to 26 pieces on 12 straight pins and 14 hairpin loops. Dimensions: 17 x 43 x 48 cm(H).
- Heathrow. **HS23243A** ..... \$105.00
- Electric, double-sided, 3 tier, Bel-Art. **188190023** ..... \$868.00

**Film & Foil**



- Parafilm**, thermoplastic sealing film
- 2" x 250' (5 cm x 76 m) Roll. **HS234526A** ..... \$21.90 ↓
  - 4" x 125' (10 cm x 38 m) Roll. **HS234526B** ..... \$22.60 ↓
  - 4" x 250' (10 cm x 76 m) Roll. **HS234526C** ..... \$39.35 ↓

Holds one or two 2" rolls, or one 4" roll.



- Parafilm Dispenser**, ABS plastic, blue. **HS234525B** ..... \$35.35



- Super-Polyfoil™**, polyester-backed aluminum foil, Bel-Art
- 4 x 4" (10.2 x 10.2 cm) Sheets, Pk/100. **132000044** ..... \$27.85
  - 5 1/8 x 5 1/8" (13 x 13 cm) Sheets, Pk/100. **132000066** ..... \$35.95