

## Filter & Blotting Papers



Ahlstrom  
Filter Paper  
USA

### Qualitative, Standard Grades (for gravity filtration)

Slow Flow (4 mL/min), 1.5 µm retention, Grade 610 (Whatman 5; S&S 602; MFS 235; Machery-Nagel 619 de; Albet FP602H; Munktell 1293)	
9 cm, Pk/100. 6100-0900 .....	\$7.95
12.5 cm, Pk/100. 6100-1250 .....	\$10.00
15 cm, Pk/100. 6100-1500 .....	\$12.10
20.5 cm, Pk/100. 6100-2050 .....	\$19.05
24 cm, Pk/100. 6100-2400 .....	\$27.05
Medium Flow (30 mL/min), 2 µm retention, Grade 642 (Whatman 2, 2CHR, 20CHR; S&S 2043B, 2045A/B; MFS 131, 232; Albet 2043B, FP593; Munktell 2, FN4, FNS)	
9 cm, Pk/100. 6420-0900 .....	\$6.20
15 cm, Pk/100. 6420-1500 .....	\$11.20
Medium Flow (55 mL/min), 2.5 µm retention, Grade 601 (Whatman 1, 1CHR; S&S 593, 597, 2043A; MFS 2; Machery-Nagel 616; Albet 2043A, FP595; Munktell 1, 1002, FN3)	
4.25 cm, Pk/100. 6010-0425 .....	\$4.75
9 cm, Pk/100. 6010-0900 .....	\$6.35
12.5 cm, Pk/100. 6010-1250 .....	\$9.25
15 cm, Pk/100. 6010-1500 .....	\$11.25
20.5 cm, Pk/100. 6010-2050 .....	\$17.60
33 cm, Pk/100. 6010-3300 .....	\$47.50
Medium Flow (60 mL/min), 6 µm retention, Grade 613 (Whatman Student; S&S 0903; Machery-Nagel 615; Munktell 3hw, 4b, 3m/N)	
9 cm, Pk/100. 6130-0900 .....	\$5.60
15 cm, Pk/100. 6130-1500 .....	\$9.80
20.5 cm, Pk/100. 6130-2050 .....	\$14.75
33 cm, Pk/100. 6130-3300 .....	\$32.75
Medium Flow (200 mL/min), 10 µm retention, Grade 631 (Whatman 4, 4CHR; S&S 604, 2040A/B; MFS 1; Albet 2040A, FP604; Munktell 1288, 1003, FN1, FN2)	
20.5 cm, Pk/100. 6310-2050 .....	\$17.75
33 cm, Pk/100. 6310-3300 .....	\$48.00
Fast Flow (235 mL/min), 25 µm retention, Grade 615 (Machery-Nagel 750 n)	
9 cm, Pk/100. 6150-0900 .....	\$5.60
12.5 cm, Pk/100. 6150-1250 .....	\$8.20
15 cm, Pk/100. 6150-1500 .....	\$9.80
20.5 cm, Pk/100. 6150-2050 .....	\$14.75
33 cm, Pk/100. 6150-3300 .....	\$32.75
Very Fast Flow (360 mL/min), 35 µm retention, Grade 617	
9 cm, Pk/50. 6170-0900 .....	\$5.10
12.5 cm, Pk/50. 6170-1250 .....	\$6.85
15 cm, Pk/50. 6170-1500 .....	\$8.00
20.5 cm, Pk/50. 6170-2050 .....	\$12.20
33 cm, Pk/50. 6170-3300 .....	\$26.10

### Qualitative, Wet-Strengthened Grades (for vacuum filtration)

Medium Flow (65 mL/min), 3 µm retention, Grade 909 (Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201; Munktell 1289, 55/N, 601/N)	
32 cm, Pk/100. 9090-3200 .....	\$51.75
Fast Flow (245 mL/min), 27 µm retention, Grade 962 (Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202; Munktell 41b, 67/N)	
11 cm, Pk/100. 9620-1100 .....	\$7.70

### Qualitative, Pleated (Fluted)

Slow Flow (30 mL/min), 2 µm retention, Grade 542 (Whatman 2V)	
12.5 cm, Pk/100. 5420-1250 .....	\$63.70
18.5 cm, Pk/100. 5420-1850 .....	\$72.10

### Quantitative, Ashless

(Designed for use in analytical and gravimetric analysis)

Slow Flow (4 mL/min), 1.5 µm retention, Grade 94 (Whatman 42; S&S 589/3; MFS 5C; Machery-Nagel 640 d; Albet FP589/3; Munktell 393)	
9 cm, Pk/100. 0940-0900 .....	\$24.95
15 cm, Pk/100. 0940-1500 .....	\$43.00
Medium Flow (30 mL/min), 2 µm retention, Grade 74 (Whatman 40; S&S 589/2, 589/4; MFS 3, 5B; Machery-Nagel 640 m; Albet 589/5; Munktell 389, 390)	
9 cm, Pk/100. 0740-0900 .....	\$24.95
15 cm, Pk/100. 0740-1500 .....	\$43.00
Fast Flow (220 mL/min), 15 µm retention, Grade 55 (Whatman 54, 541; S&S 1505; Albet FP1505, FP1573; Munktell 388)	
9 cm, Pk/100. 0550-0900 .....	\$28.20
15 cm, Pk/100. 0550-1500 .....	\$45.30

### Blotting/Chromatography/Electrophoresis Papers

0.41 mm thick, 150 mm/30 min absorption, Grade 237 (Whatman 3, 3CHR; S&S 22/2, 591, 2469, 5703; Albet 2040B, FP591; Munktell 3S/H, FN7a)	
46 x 57 cm, Pk/100. 2378-4657 .....	\$198.00 ↓
1.02 mm thick, 180 mm/30 min absorption, Grade 243 (Whatman 17CHR; S&S 2668; Munktell 151)	
11 x 14 cm, Pk/100. 2438-1114 .....	\$17.15
20 x 20 cm, Pk/100. 2438-2020 .....	\$33.55
46 x 57 cm, Pk/100. 2438-4657 .....	\$289.00 ↓
2.45 mm thick, 240 mm/30 min absorption, Grade 320 (S&S 2727)	
20 x 20 cm, Pk/50. 3208-2020 .....	\$54.70



Smooth tips &  
polypropylene grips

Filter Forceps, 11.5 cm. 300-025 .....	\$5.80
--	--------