

Hot Plate/Magnetic Stirrer/Stirring Hot Plate



Hot Plate



Magnetic Stirrer



Stirring Hot Plate

- Hot Plate**, stainless steel, 50-300°C. **DB-1** \$173.00
- Magnetic Stirrer**, 7.5 x 7.5" platform, 60-1500 rpm. **H4000-S** ... \$199.00
- Stirring Hot Plate**, external temperature probe. **88-1** \$272.00



Ice Buckets




Ice Pans

Ice Buckets with lids, expanded PVC foam, Bel-Art, **USA**

- 2.5 L Capacity
 - Black. **188482002** \$63.50
 - Blue. **188482001** \$63.50
 - Purple. **188482006** \$63.50
- 4 L Capacity
 - Black. **188484002** \$75.50
 - Blue. **188484001** \$75.50
 - Purple. **188484006** \$75.50

Ice Pans, expanded PVC foam, Bel-Art, **USA**

- 1 L Capacity, 15.5 x 15.5 x 8.5 cm(H)
 - Black. **188481102** \$44.65
 - Blue. **188481101** \$44.65
 - Purple. **188481106** \$44.65
- 4 L Capacity, 33 x 24.5 x 12.5 cm(H)
 - Black. **188484102** \$75.50
 - Blue. **188484101** \$75.50
 - Purple. **188484106** \$75.50
- 9 L Capacity, 44 x 35 x 13 cm(H)
 - Black. **188489102** \$90.95
 - Blue. **188489101** \$90.95
 - Purple. **188489106** \$90.95




Tube Racks, Floating, PVC Foam, Bel-Art, **USA**

- Circular. *Recommended for ice buckets*
 - 44-Hole, 1.5 mL tubes, 22 cm dia. x 1 cm \$28.05
Black/**188480012** Green/**188480014** Red/**188480013**
 - 30-Hole, 14/15 mL tubes, 22 cm dia. x 1.7 cm \$29.95
Black/**188480022** Green/**188480024** Red/**188480023**
- Rectangular. *Recommended for ice pans*
 - 36-Hole, 1.5 mL tubes, 23.5 x 11 x 1 cm \$28.05
Black/**188480112** Green/**188480114** Red/**188480113**
 - 24-Hole, 14/15 mL tubes, 23.5 x 11 x 1.7 cm \$29.95
Black/**188480122** Green/**188480124** Red/**188480123**




Replace aluminum blocks, water, and ice in baths and buckets with this thermal metallic media.



Lab Armor™ Beads, USA

- 0.75 L (replaces two blocks in a dry bath). **42370-750** \$98.00
- 2 L. **42370-002** \$191.00
(B-0700-2L)
- 4 L. **42370-004** \$314.00
(B-0700-4L)
- 8 L (2 x 4 L). **42370-008** \$490.00
(B-0700-8L)

Label plastics, glass and metals with these chemically inert and temperature resistant (-196 to 150°C) polyolefin tags.



Labels


- Cryo-Babies®, white, Diversified Biotech
 - 24 x 13 mm (0.94 x 0.5"), Roll/1000. **CRY-0500** \$29.75
Fits 0.5 mL tubes.
 - 33 x 13 mm (1.28 x 0.5"), Roll/1000. **CRY-1000** \$29.75
Fits 1.5-2.0 mL tubes.

Kraft paper with rubber resin adhesive.



Label Tape, Write-On™

- 1" Core, 15 yard roll
 - ½" Wide \$5.00
 - Green. **134623005** Yellow. **134622005**
 - White. **134620005**
 - ¾" Wide \$5.70
 - Green. **134623075** Yellow. **134622075**
 - White. **134620075**
- 3" Core, 40 yard roll
 - ½" Wide \$9.70
 - Green. **134633005** White. **134630005**
 - Red. **134631005** Yellow. **134632005**
 - ¾" Wide \$12.85
 - Green. **134633075** White. **134630075**
 - Red. **134631075** Yellow. **134632075**



- Label Tape Dispenser**, economy, Bel-Art. **134730000** \$33.70

Tubes



- Centrifuge, conical, polypropylene, sterile
15 mL, sterile bag/50, Cs/500. **HS4426B** \$83.20
- 50 mL, max RCF is 9,500 x g
Natural, sterile bag/50, Cs/500. **HS4427B** \$114.00
- Black, sterile bag/50, Cs/500. **HS4429** \$165.00



- Centrifuge w/Attached Cap, polypropylene, max RCF is 5,000 x g
5 mL, 0.5 mL graduations, green, Pk/250. **T2079** \$75.95



Polystyrene with dual-position caps.

- Culture, 12 x 75 mm (5 mL), sterile, Simport
Cs/500, sterile bag/25. **T415-2** \$63.00
- Cs/1000, sterile bag/125. **T415-3** \$94.00



Natural



Black

- Microcentrifuge, Economy, polypropylene (PP), natural color
0.5 mL, 0.1 mL graduations, natural, Pk/1000. **HS4422** \$32.20
- 1.5 mL, 0.25 mL graduations
Natural, max RCF is 17,000 x g
Pk/500. **HS4323** \$15.05
- Pk/1000. **HS4325** \$29.95
- Black, max RCF is 18,000 x g, Pk/500. **HS4323K** \$17.40



0.6 mL

1.5 mL

2.0 mL

Tubes will withstand centrifugal forces over 30,000 x g and temperatures from -80° to 121°C. Flat, textured caps are pierceable.

- Microcentrifuge, Premium, PP, G-Tube™, Pk/500, USA
Amber color protects light-sensitive specimens.
- 0.6 mL \$24.65
- Amber. **4048** Green. **4044** Natural. **4040**
- Blue. **4043** Lavender. **4047** Yellow. **4042**
- Mixed (blue, green, orange, pink, & yellow). **4049** \$29.70
- 1.5 mL, 0.5 mL graduations, frosted area on side \$24.65
- Amber. **4038** Green. **4034** Natural. **4030**
- Blue. **4033** Lavender. **4037** Yellow. **4032**
- Mixed (blue, green, orange, pink, & yellow). **4039** \$29.70

Continued

- 2.0 mL, 0.5 mL graduations, frosted area on side \$24.10
- Amber. **4058** Green. **4054** Natural. **4050**
- Blue. **4053** Lavender. **4057** Yellow. **4052**
- Mixed (blue, green, orange, pink, & yellow). **4059** \$31.55

Microcentrifuge, Siliconized, G-Tube™, Pk/500, USA

- 1.5 mL, natural color. **4165SL** \$46.00
(Siliconized products should NOT be autoclaved)



Securely close lids on 0.5 and 1.5 mL microcentrifuge tubes with these autoclavable, polypropylene holders. Extract tubes from centrifuges, dry baths, racks, etc. using large tabs. Assorted colors.


- Microcentrifuge Tube Holders, Pk/100. **HS4350** \$21.35

PCR, 0.2 mL thin-wall, polypropylene, USA

Strips (*caps are sold separately; available in natural only*)

- 8-Tube Strips, Pk/125 \$82.25
- Blue. **5010-4** Natural. **5010-1**
- 8-Tube Cap Strips, domed, Pk/300. **5030-1** \$80.55
- 12-Tube Strips, Pk/100 \$100.00
- Blue. **5020-4** Natural. **5020-1**
- 12-Tube Cap Strips, domed, Pk/200. **5040-1** \$80.55

Tube Racks



Floating, PVC Foam, Bel-Art, USA

Circular. *Recommended for ice buckets*

- 44-Hole, 1.5 mL tubes, 22 cm dia. x 1 cm \$28.05
- Black/**188480012** Green/**188480014** Red/**188480013**
- 30-Hole, 14/15 mL tubes, 22 cm dia. x 1.7 cm \$29.95
- Black/**188480022** Green/**188480024** Red/**188480023**

Rectangular. *Recommended for ice pans*


- 36-Hole, 1.5 mL tubes, 23.5 x 11 x 1 cm \$28.05
- Black/**188480112** Green/**188480114** Red/**188480113**
- 24-Hole, 14/15 mL tubes, 23.5 x 11 x 1.7 cm \$29.95
- Black/**188480122** Green/**188480124** Red/**188480123**



Floating, HDPE Foam w/Handle, Pk/5

- 18-Hole, 0.2/0.5/1.5-2.0 mL tubes, blue circle. **HS2166** \$13.60
- 24-Hole, 1.5-2.0 mL tubes, yellow rectangle. **HS2165C** \$13.40
- 8-Hole, 15 mL tubes, blue parallelogram. **HS2165B** \$11.80
- 4-Hole, 50 mL tubes, green diamond. **HS2165A** \$13.20

Continued



Floating, Polyester Foam, gray
30-Hole, 250/400 µL tubes, Pk/6. **HS21645C**..... \$20.50




Grip racks, polypropylene, No-Wire™, Bel-Art
90-Place, 10-13 mm tubes, each. **187490000** \$21.10
60-Place, 15-16 mm tubes, each. **187490001** \$21.10
40-Place, 18-20 mm tubes, each. **187490002** \$21.10



Cube Rack



4-Way Rack



2-Tier Rack

4-Way Rotatable Rack

Interlocking racks, polypropylene, mixed colors
Cube, Pk/5. **HS29050A** \$36.55
Holds (4) 50 mL tubes, (10) 15 mL tubes,
(12) 12 x 75 or 12 x 100 tubes, or (16) 1.5-2.0 mL tubes.
2-Tier, red, Pk/6. **HS4284B** \$40.40
Holds (20) 0.5-2.0 mL tubes on two levels.
4-Way, Pk/5. **HS29022G** \$35.65
Holds (4) 50 mL tubes, (12) 15 mL tubes,
(32) 1.5 mL tubes, or (32) 0.5 mL tubes.
4-Way Rotatable, Rota-Rack™, four connected modules
Large Tubes, 22.4 x 11 x 11 cm. **HS29040A** \$36.90
Each module holds (2) 50 mL tubes, (8) 15 mL tubes,
(6) 20 mm tubes, or (10) 10-13 mm tubes.
Microtubes, 21 x 10 x 10 cm. **HS29040B** \$36.90
Each module holds (6) 15 mL tubes, (9) 1.5/2.0 mL tubes,
(12) 0.5/0.6 mL tubes, or (32) 0.2 mL PCR tubes or
(4) 8-tube strips.

HS2860A/H



Microtube Storage Boxes

Cardboard, vented
2" Box Only, 5¼ x 5¼ x 2", Pk/12. **HS2860A** \$24.30
2" box w/dividers ea pk/12
100-Cell (10 x 10). **HS2860A/E** \$4.05 44.20
81-Cell (9 x 9). **HS2860A/H** \$4.05 44.20
64-Cell (8 x 8). **HS2860A/F** \$4.00 43.60
2" Mini box w/dividers, 3 x 3 x 2"
25-Cell (5 x 5), Pk/12. **MB-2/MD-25** \$45.60
Holds 13 mm vials and 4.0-5.0 mL cryovials.
3" Box Only, 5¼ x 5¼ x 3", Pk/12. **HS2860B** \$26.30
3" box w/dividers ea pk/12
100-Cell (10 x 10). **HS2860B/E** \$4.15 45.40
81-Cell (9 x 9). **HS2860B/H** \$4.15 45.40
64-Cell (8 x 8). **HS2860B/F** \$4.10 44.80
Dividers only, fit standard 2 & 3" boxes
100-Cell (10 x 10), Pk/12. **HS2860E** \$23.20
Holds 1.5-2.0 mL tubes and 1.2-2.0 mL cryovials.
81-Cell (9 x 9), Pk/12. **HS2860H** \$23.20
Holds 13 mm vials and 4.0-5.0 mL cryovials.
64-Cell (8 x 8), Pk/12. **HS2860F** \$22.05
Holds 14 mm vials.
25-Cell (5 x 5), 25 mm cell, Pk/12. **D-25**..... \$21.60
16-Cell (4 x 4), 31 mm cell, Pk/12. **D-16**..... \$21.60



HS23381B



188490002

Polycarbonate, store vials from -196°C to +121°C
Arctic Squares®
25-Well (5 x 5), friction-fit lid, Pk/8. **HS23325A** \$49.65
Holds vials up to 2 mL; red. 7.6 x 7.6 x 5.3 cm
81-Well (9 x 9), hinged lid, Pk/4. **HS23381B** \$40.60
Holds vials up to 2 mL; mixed. 13.3 x 13.3 x 5.3 cm
Cryo-Safe™, 13.3 x 13.3 x 9.5 cm, Bel-Art
81-Well (9 x 9), Pk/4. **188490002** \$48.95
Holds 3.0/5.0 mL vials.




HS23255



HS1065

Polypropylene
81-Well (9 x 9), 13 x 13 x 4.7 cm
Holds 1.5-2.0 mL tubes and cryovials.
Features imprinted coordinates on a colored base and
a translucent, friction-fit lid with frosted panels.
Blue Base, each. **HS23255B**..... \$6.40
Pink Base, each. **HS23255D** \$6.40
100-Well (10 x 10), assorted colors, Pk/5. **HS1065**..... \$31.40
Holds 1.5-2.0 mL tubes and cryovials. 14.1 x 15.1 x 5.7 cm



HS29025



HS2345

Microtube Racks

Polypropylene

- 80-Well (16 x 5), molded-in coordinates, 22.5 x 6.7 x 2.8 cm
Holds 1.5-2.0 mL tubes in conical wells.
- Blue, Pk/5. **HS29025B** \$23.75
- Pink, Pk/5. **HS29025D** \$23.75
- 96-Well reversible (12 x 8), snap-on cover, 24.6 x 12.1 x 5 cm
Holds 1.5-2.0 mL tubes on one side or 0.5 mL on reverse.
- Blue, each. **HS2345C** \$13.50
- Green, each. **HS2345D** \$13.50
- Pink, each. **HS2345E** \$13.50


Holds 1.5-2.0 mL tubes
in a two-piece rack.



- Polystyrene Foam (styrofoam)
- 50-Well, Pk/10. **HS4283P** \$36.90
- 100-Well, Pk/10. **HS4282** \$46.45



- PCR Reversible racks, assorted colors, Pk/5. **HS2344A** \$47.95
- Holds 0.2, 0.5 and 1.5 mL tubes.



- Peg racks, polypropylene
- 96-Place, 10-13 mm tubes, blue, each. **HS24312B** \$11.80
- 50-Place, 14-17 mm tubes, blue, each. **HS24311B** \$10.70



- Refrigerator Racks, two-piece polystyrene
- 50-Well, 1.5 mL tubes, each. **HS28230A** \$35.35
- 15-Well, 15/50 mL tubes, each. **HS28230B** \$35.35



- Slant racks, 40-place, polypropylene, Bel-Art
- Holds 13-16 mm Tubes, each. **188900016** \$46.75
- Holds 16-20 mm Tubes, each. **188900020** \$46.75




Tilt your test tube rack 25°, 35°, or 45°. Create agar slant tubes, or adjust the angle of tube racks when pipetting to reduce hand strain.

- Tube Rack Slanter. **188900100** \$15.95
- Platform size: 19.6 cm (W) x 12.7 cm (D)



- Stackable, 72-place, polypropylene, Bel-Art
- Holds 1.5 mL Tubes, each. **188410000** \$32.80



- SmoothRack™, 72-place, polypropylene, blue, USA
- Holds 10-17 mm Tubes, each. **SR-B** \$14.75



Remains stable
in water bath.

- Snap-Together racks, polypropylene, blue
- 90-Place, 13 mm tubes, each. **HS24305A** \$12.70
- 60-Place, 16 mm tubes, each. **HS24305B** \$12.70
- 40-Place, 19 mm tubes, each. **HS24305C** \$12.70