

—Hot Plate/Magnetic Stirrer/Stirring Hot Plate—



Hot Plate



Magnetic Stirrer




Stirring Hot Plate

- Hot Plate**, stainless steel, 50-300°C. **DB-1** ..... \$173.00
- Magnetic Stirrer**, 7.5 x 7.5" platform, 60-1500 rpm. **H4000-S** ... \$199.00
- Stirring Hot Plate**, external temperature probe. **88-1** ..... \$272.00



Ice Buckets




Ice Pans

**Ice Buckets** with lids, expanded PVC foam, Bel-Art, **USA**

- 2.5 L Capacity
  - Black. **188482002** ..... \$63.50
  - Blue. **188482001** ..... \$63.50
  - Purple. **188482006** ..... \$63.50
- 4 L Capacity
  - Black. **188484002** ..... \$75.50
  - Blue. **188484001** ..... \$75.50
  - Purple. **188484006** ..... \$75.50

**Ice Pans**, expanded PVC foam, Bel-Art, **USA**

- 1 L Capacity, 15.5 x 15.5 x 8.5 cm(H)
  - Black. **188481102** ..... \$44.65
  - Blue. **188481101** ..... \$44.65
  - Purple. **188481106** ..... \$44.65
- 4 L Capacity, 33 x 24.5 x 12.5 cm(H)
  - Black. **188484102** ..... \$75.50
  - Blue. **188484101** ..... \$75.50
  - Purple. **188484106** ..... \$75.50
- 9 L Capacity, 44 x 35 x 13 cm(H)
  - Black. **188489102** ..... \$90.95
  - Blue. **188489101** ..... \$90.95
  - Purple. **188489106** ..... \$90.95




**Tube Racks**, Floating, PVC Foam, Bel-Art, **USA**

- Circular. *Recommended for ice buckets*
  - 44-Hole, 1.5 mL tubes, 22 cm dia. x 1 cm ..... \$28.05  
Black/**188480012** Green/**188480014** Red/**188480013**
  - 30-Hole, 14/15 mL tubes, 22 cm dia. x 1.7 cm ..... \$29.95  
Black/**188480022** Green/**188480024** Red/**188480023**
- Rectangular. *Recommended for ice pans*
  - 36-Hole, 1.5 mL tubes, 23.5 x 11 x 1 cm ..... \$28.05  
Black/**188480112** Green/**188480114** Red/**188480113**
  - 24-Hole, 14/15 mL tubes, 23.5 x 11 x 1.7 cm ..... \$29.95  
Black/**188480122** Green/**188480124** Red/**188480123**




Replace aluminum blocks, water, and ice in baths and buckets with this thermal metallic media.



**Lab Armor™ Beads, USA**

- 0.75 L (replaces two blocks in a dry bath). **42370-750** ..... \$98.00
- 2 L. **42370-002** ..... \$191.00  
(B-0700-2L)
- 4 L. **42370-004** ..... \$314.00  
(B-0700-4L)
- 8 L (2 x 4 L). **42370-008** ..... \$490.00  
(B-0700-8L)

Label plastics, glass and metals with these chemically inert and temperature resistant (-196 to 150°C) polyolefin tags.



**Labels**


- Cryo-Babies®, white, Diversified Biotech
  - 24 x 13 mm (0.94 x 0.5"), Roll/1000. **CRY-0500** ..... \$29.75  
Fits 0.5 mL tubes.
  - 33 x 13 mm (1.28 x 0.5"), Roll/1000. **CRY-1000** ..... \$29.75  
Fits 1.5-2.0 mL tubes.

Kraft paper with rubber resin adhesive.



**Label Tape**, Write-On™

- 1" Core, 15 yard roll
  - ½" Wide ..... \$5.00
    - Green. **134623005**      Yellow. **134622005**
    - White. **134620005**
  - ¾" Wide ..... \$5.70
    - Green. **134623075**      Yellow. **134622075**
    - White. **134620075**
- 3" Core, 40 yard roll
  - ½" Wide ..... \$9.70
    - Green. **134633005**      White. **134630005**
    - Red. **134631005**      Yellow. **134632005**
  - ¾" Wide ..... \$12.85
    - Green. **134633075**      White. **134630075**
    - Red. **134631075**      Yellow. **134632075**



- Label Tape Dispenser**, economy, Bel-Art. **134730000** ..... \$33.70

**Tubes**



- Centrifuge, conical, polypropylene, sterile  
15 mL, sterile bag/50, Cs/500. **HS4426B** ..... \$83.20
- 50 mL, max RCF is 9,500 x g  
Natural, sterile bag/50, Cs/500. **HS4427B** ..... \$114.00
- Black, sterile bag/50, Cs/500. **HS4429** ..... \$165.00



- Centrifuge w/Attached Cap, polypropylene, max RCF is 5,000 x g  
5 mL, 0.5 mL graduations, green, Pk/250. **T2079** ..... \$75.95



Polystyrene with dual-position caps.

- Culture, 12 x 75 mm (5 mL), sterile, Simport  
Cs/500, sterile bag/25. **T415-2** ..... \$63.00
- Cs/1000, sterile bag/125. **T415-3** ..... \$94.00




Natural



Black

- Microcentrifuge, Economy, polypropylene (PP), natural color  
0.5 mL, 0.1 mL graduations, natural, Pk/1000. **HS4422** ..... \$32.20
- 1.5 mL, 0.25 mL graduations  
Natural, max RCF is 17,000 x g  
Pk/500. **HS4323** ..... \$15.05
- Pk/1000. **HS4325** ..... \$29.95
- Black, max RCF is 18,000 x g, Pk/500. **HS4323K** ..... \$17.40



0.6 mL

1.5 mL

2.0 mL

Tubes will withstand centrifugal forces over 30,000 x g and temperatures from -80° to 121°C. Flat, textured caps are pierceable.

- Microcentrifuge, Premium, PP, G-Tube™, Pk/500, USA  
*Amber color protects light-sensitive specimens.*
- 0.6 mL ..... \$24.65
- Amber. **4048**    Green. **4044**    Natural. **4040**
- Blue. **4043**    Lavender. **4047**    Yellow. **4042**
- Mixed (blue, green, orange, pink, & yellow). **4049** ..... \$29.70
- 1.5 mL, 0.5 mL graduations, frosted area on side ..... \$24.65
- Amber. **4038**    Green. **4034**    Natural. **4030**
- Blue. **4033**    Lavender. **4037**    Yellow. **4032**
- Mixed (blue, green, orange, pink, & yellow). **4039** ..... \$29.70

Continued

- 2.0 mL, 0.5 mL graduations, frosted area on side ..... \$24.10
- Amber. **4058**    Green. **4054**    Natural. **4050**
- Blue. **4053**    Lavender. **4057**    Yellow. **4052**
- Mixed (blue, green, orange, pink, & yellow). **4059** ..... \$31.55

- Microcentrifuge, Siliconized, G-Tube™, Pk/500, USA  
1.5 mL, natural color. **4165SL** ..... \$46.00  
*(Siliconized products should NOT be autoclaved)*




Securely close lids on 0.5 and 1.5 mL microcentrifuge tubes with these autoclavable, polypropylene holders. Extract tubes from centrifuges, dry baths, racks, etc. using large tabs. Assorted colors.

- Microcentrifuge Tube Holders, Pk/100. **HS4350** ..... \$21.35

- PCR, 0.2 mL thin-wall, polypropylene, USA  
Strips (*caps are sold separately; available in natural only*)
- 8-Tube Strips, Pk/125 ..... \$82.25
- Blue. **5010-4**    Natural. **5010-1**
- 8-Tube Cap Strips, domed, Pk/300. **5030-1** ..... \$80.55
- 12-Tube Strips, Pk/100 ..... \$100.00
- Blue. **5020-4**    Natural. **5020-1**
- 12-Tube Cap Strips, domed, Pk/200. **5040-1** ..... \$80.55

**Tube Racks**




- Floating, PVC Foam, Bel-Art, USA
- Circular. *Recommended for ice buckets*
- 44-Hole, 1.5 mL tubes, 22 cm dia. x 1 cm ..... \$28.05
- Black/**188480012**    Green/**188480014**    Red/**188480013**
- 30-Hole, 14/15 mL tubes, 22 cm dia. x 1.7 cm ..... \$29.95
- Black/**188480022**    Green/**188480024**    Red/**188480023**
- Rectangular. *Recommended for ice pans*
- 36-Hole, 1.5 mL tubes, 23.5 x 11 x 1 cm ..... \$28.05
- Black/**188480112**    Green/**188480114**    Red/**188480113**
- 24-Hole, 14/15 mL tubes, 23.5 x 11 x 1.7 cm ..... \$29.95
- Black/**188480122**    Green/**188480124**    Red/**188480123**



- Floating, HDPE Foam w/Handle, Pk/5
- 18-Hole, 0.2/0.5/1.5-2.0 mL tubes, blue circle. **HS2166** ..... \$13.60
- 24-Hole, 1.5-2.0 mL tubes, yellow rectangle. **HS2165C** ..... \$13.40
- 8-Hole, 15 mL tubes, blue parallelogram. **HS2165B** ..... \$11.80
- 4-Hole, 50 mL tubes, green diamond. **HS2165A** ..... \$13.20

Continued



**Floating, Polyester Foam, gray**  
30-Hole, 250/400 µL tubes, Pk/6. **HS21645C**..... \$20.50




**Grip racks, polypropylene, No-Wire™, Bel-Art**  
90-Place, 10-13 mm tubes, each. **187490000** ..... \$21.10  
60-Place, 15-16 mm tubes, each. **187490001** ..... \$21.10  
40-Place, 18-20 mm tubes, each. **187490002** ..... \$21.10



Cube Rack



4-Way Rack



2-Tier Rack

4-Way Rotatable Rack

**Interlocking racks, polypropylene, mixed colors**  
Cube, Pk/5. **HS29050A** ..... \$36.55  
Holds (4) 50 mL tubes, (10) 15 mL tubes,  
(12) 12 x 75 or 12 x 100 tubes, or (16) 1.5-2.0 mL tubes.  
2-Tier, red, Pk/6. **HS4284B** ..... \$40.40  
Holds (20) 0.5-2.0 mL tubes on two levels.  
4-Way, Pk/5. **HS29022G** ..... \$35.65  
Holds (4) 50 mL tubes, (12) 15 mL tubes,  
(32) 1.5 mL tubes, or (32) 0.5 mL tubes.  
4-Way Rotatable, Rota-Rack™, four connected modules  
Large Tubes, 22.4 x 11 x 11 cm. **HS29040A** ..... \$36.90  
Each module holds (2) 50 mL tubes, (8) 15 mL tubes,  
(6) 20 mm tubes, or (10) 10-13 mm tubes.  
Microtubes, 21 x 10 x 10 cm. **HS29040B** ..... \$36.90  
Each module holds (6) 15 mL tubes, (9) 1.5/2.0 mL tubes,  
(12) 0.5/0.6 mL tubes, or (32) 0.2 mL PCR tubes or  
(4) 8-tube strips.

HS2860A/H



**Microtube Storage Boxes**

Cardboard, vented  
2" Box Only, 5¼ x 5¼ x 2", Pk/12. **HS2860A** ..... \$24.30  
2" box w/dividers ea pk/12  
100-Cell (10 x 10). **HS2860A/E** ..... \$4.05 44.20  
81-Cell (9 x 9). **HS2860A/H** ..... \$4.05 44.20  
64-Cell (8 x 8). **HS2860A/F** ..... \$4.00 43.60  
2" Mini box w/dividers, 3 x 3 x 2"  
25-Cell (5 x 5), Pk/12. **MB-2/MD-25** ..... \$45.60  
Holds 13 mm vials and 4.0-5.0 mL cryovials.  
3" Box Only, 5¼ x 5¼ x 3", Pk/12. **HS2860B** ..... \$26.30  
3" box w/dividers ea pk/12  
100-Cell (10 x 10). **HS2860B/E** ..... \$4.15 45.40  
81-Cell (9 x 9). **HS2860B/H** ..... \$4.15 45.40  
64-Cell (8 x 8). **HS2860B/F** ..... \$4.10 44.80  
Dividers only, fit standard 2 & 3" boxes  
100-Cell (10 x 10), Pk/12. **HS2860E** ..... \$23.20  
Holds 1.5-2.0 mL tubes and 1.2-2.0 mL cryovials.  
81-Cell (9 x 9), Pk/12. **HS2860H** ..... \$23.20  
Holds 13 mm vials and 4.0-5.0 mL cryovials.  
64-Cell (8 x 8), Pk/12. **HS2860F** ..... \$22.05  
Holds 14 mm vials.  
25-Cell (5 x 5), 25 mm cell, Pk/12. **D-25**..... \$21.60  
16-Cell (4 x 4), 31 mm cell, Pk/12. **D-16**..... \$21.60



HS23381B



188490002

Polycarbonate, store vials from -196°C to +121°C  
Arctic Squares®  
25-Well (5 x 5), friction-fit lid, Pk/8. **HS23325A** ..... \$49.65  
Holds vials up to 2 mL; red. 7.6 x 7.6 x 5.3 cm  
81-Well (9 x 9), hinged lid, Pk/4. **HS23381B** ..... \$40.60  
Holds vials up to 2 mL; mixed. 13.3 x 13.3 x 5.3 cm  
Cryo-Safe™, 13.3 x 13.3 x 9.5 cm, Bel-Art  
81-Well (9 x 9), Pk/4. **188490002** ..... \$48.95  
Holds 3.0/5.0 mL vials.




HS23255



HS1065

Polypropylene  
81-Well (9 x 9), 13 x 13 x 4.7 cm  
Holds 1.5-2.0 mL tubes and cryovials.  
Features imprinted coordinates on a colored base and  
a translucent, friction-fit lid with frosted panels.  
Blue Base, each. **HS23255B**..... \$6.40  
Pink Base, each. **HS23255D** ..... \$6.40  
100-Well (10 x 10), assorted colors, Pk/5. **HS1065**..... \$31.40  
Holds 1.5-2.0 mL tubes and cryovials. 14.1 x 15.1 x 5.7 cm



HS29025



HS2345

**Microtube Racks**

Polypropylene

- 80-Well (16 x 5), molded-in coordinates, 22.5 x 6.7 x 2.8 cm  
Holds 1.5-2.0 mL tubes in conical wells.
- Blue, Pk/5. **HS29025B** ..... \$23.75
- Pink, Pk/5. **HS29025D** ..... \$23.75
- 96-Well reversible (12 x 8), snap-on cover, 24.6 x 12.1 x 5 cm  
Holds 1.5-2.0 mL tubes on one side or 0.5 mL on reverse.
- Blue, each. **HS2345C** ..... \$13.50
- Green, each. **HS2345D** ..... \$13.50
- Pink, each. **HS2345E** ..... \$13.50

Holds 1.5-2.0 mL tubes  
in a two-piece rack.



Polystyrene Foam (styrofoam)

- 50-Well, Pk/10. **HS4283P** ..... \$36.90
- 100-Well, Pk/10. **HS4282** ..... \$46.45



- PCR Reversible racks**, assorted colors, Pk/5. **HS2344A** ..... \$47.95
- Holds 0.2, 0.5 and 1.5 mL tubes.



**Peg racks**, polypropylene

- 96-Place, 10-13 mm tubes, blue, each. **HS24312B** ..... \$11.80
- 50-Place, 14-17 mm tubes, blue, each. **HS24311B** ..... \$10.70



**Refrigerator Racks**, two-piece polystyrene

- 50-Well, 1.5 mL tubes, each. **HS28230A** ..... \$35.35
- 15-Well, 15/50 mL tubes, each. **HS28230B** ..... \$35.35



**Slant racks**, 40-place, polypropylene, Bel-Art

- Holds 13-16 mm Tubes, each. **188900016** ..... \$46.75
- Holds 16-20 mm Tubes, each. **188900020** ..... \$46.75




Tilt your test tube rack 25°, 35°, or 45°. Create agar slant tubes, or adjust the angle of tube racks when pipetting to reduce hand strain.

- Tube Rack Slanter**. **188900100** ..... \$15.95
- Platform size: 19.6 cm (W) x 12.7 cm (D)



**Stackable**, 72-place, polypropylene, Bel-Art

- Holds 1.5 mL Tubes, each. **188410000** ..... \$32.80



**SmoothRack™**, 72-place, polypropylene, blue, USA


- Holds 10-17 mm Tubes, each. **SR-B** ..... \$14.75



Remains stable  
in water bath.

**Snap-Together racks**, polypropylene, blue

- 90-Place, 13 mm tubes, each. **HS24305A** ..... \$12.70
- 60-Place, 16 mm tubes, each. **HS24305B** ..... \$12.70
- 40-Place, 19 mm tubes, each. **HS24305C** ..... \$12.70




**Teardrop racks, OneRack™, POM**

**Half-Size**


- 36-Place, holds 16 mm tubes, 12.7 x 12.7 x 7 cm(H)
- Blue, Pk/8. **HS27502B** ..... \$92.00
- Red, Pk/8. **HS27502D** ..... \$92.00

**Full-Size**

- 72-Place, holds 13 mm tubes, 20.2 x 10.4 x 5.9 cm(H)
- Blue, Pk/8. **HS27511B** ..... \$120.00
- Red, Pk/8. **HS27511D** ..... \$120.00



- 8-Place, epoxy-coated wire, each. **HS232461** ..... \$24.40
- 9-Place (3 x 3), stainless steel, USA, each. **50ML-9** ..... \$14.80
- 16-Place, epoxy-coated wire, each. **187940001** ..... \$38.60
- 30-Place (3 x 10), stainless steel, USA, each. **50ML-30** ..... \$33.40
- Dimensions: 40.7 x 12.7 x 10.2 cm(H).
- 42-Place (3 x 14), stainless steel, USA, each. **50ML-42** ..... \$43.00
- Dimensions: 54.6 x 12.7 x 10.2 cm(H).




120-Place Wire Rack

**Wire racks, epoxy-coated, blue, Poxygrip®, Bel-Art**

- 120-Place (6 x 20), 10-13 mm tubes, ea. **187860750** ..... \$34.35
- Dimensions: 33.1 x 11 x 6 cm(H).
- 100-Place (5 x 20), 15-16 mm tubes, ea. **187860780** ..... \$34.35
- Dimensions: 38.6 x 10.8 x 6.4 cm(H).
- 180-Place (9 x 20), 15-16 mm tubes, ea. **187860770** ..... \$35.35
- Dimensions: 38.6 x 18.3 x 6.4 cm(H).

HS23215



15ML-25

**15 mL tube racks**


- 15-Place, epoxy-coated wire, each. **HS23215** ..... \$24.15
- 25-Place (5 x 5), stainless steel, USA, each. **15ML-25** ..... \$14.80
- 80-Place (5 x 16), stainless steel, USA, each. **15ML-80** ..... \$33.40
- Dimensions: 40.7 x 12.7 x 10.2 cm(H).
- 110-Place (5 x 22), stainless steel, USA, each. **15ML-110** ..... \$43.00
- Dimensions: 54.6 x 12.7 x 10.2 cm(H).



Wire



Friction-Fit



Interlocking


**50 mL tube racks**

- 1-Place
- Wire, epoxy-coated, each. **HS23051** ..... \$10.45
- Friction-Fit, polystyrene, Pk/5. **HS23053A** ..... \$16.10
- Interlocking, polypropylene, Pk/5. **HS23052** ..... \$31.55



EV432

Polypropylene vacuum bottles are autoclavable and withstand 700 mm Hg systems.



199170001

**Vacuum Waste Collection Systems**

- 1.9 L (0.5 gal), Bel-Art. **199170001** ..... \$165.00
- Attached trap protects vacuum source from overflow.
- 4 L (1 gal), Argos. **EV432** ..... \$122.00
- Cap has traditional barbed fitting.



Transport and store biological specimens in these self-standing tubes. External threads minimize contamination. Dimensions of tube with cap: 17 mm diameter x 84 mm high.

**Vials, polypropylene with cap, 10 mL, Simport**

- Graduated, Cs/1000. **T550-10ATPR/T550W** ..... \$151.00
- Tubes have a writing surface and 0.5 mL graduations from 1-9 mL.
- Plain, Cs/1000. **T550-10AT/T550W** ..... \$120.00



**Vial Decapping Pliers, vinyl grip, USA**

- For 13 mm Crimp Seals. **D-13** ..... \$55.25 ↓
- For 20 mm Crimp Seals. **D-20** ..... \$55.25 ↓



Square



Syracuse

**Watch Glasses, flint glass, USA**

- Square, ground front, 41 mm. **910-B** ..... \$7.80
- Syracuse, grooved for stacking, 67 mm OD. **913-A** ..... \$5.75