



Pipet Fillers

- 3-Valve, rubber, black. **207A** \$14.40
- Pipette Pump™, wheel-type, Bel-Art
 - Fits Pipets Up To 2 mL, blue. **378970000** \$21.85
 - Fits Pipets Up To 10 mL, green. **378980000** \$22.85
 - Fits Pipets Up To 25 mL, red. **378990000** \$31.55
- Pipette Pump™ II, wheel-type, quick-release ea pk/12
 - Fits Pipets Up To 2 mL, blue. **379111002** \$23.40 253.00
 - Fits Pipets Up To 10 mL, green. **379111010** \$23.65 256.00
 - Fits Pipets Up To 25 mL, red. **379111025** \$31.55 340.00




Motorized

- Ergopet® I Memory, Bel-Art. **379100011** \$632.00
- Replacement Filters, Pk/12. **379080251** \$93.45



Dimensions: 29 x 8.8 x 40 cm

- Pipet Rack, magnetized, blue. HS20612M** \$66.25




Store up to 4 sterilization cans with this epoxy-coated wire rack. (Pipet cans not included)

- Pipet Sterilizing Box Rack, Bel-Art. 189420000** \$37.90



- Pipet Tip Disposal Box. HS234634** \$80.00



Pipets, Transfer, polyethylene Nonsterile

- 5 mL Capacity (1 mL graduated), 25 drops/mL
 - Bg/25. **PTP-02-025** \$1.00
 - Bx/500. **PTP-02-500** \$15.55
 - Cs/10 Bx. **PTP-02-5000** \$140.00
- 7 mL Capacity (3 mL graduated), 14 drops/mL
 - Bg/25. **PTP-01-025** \$1.00
 - Bx/500. **PTP-01-500** \$15.55
 - Cs/10 Bx. **PTP-01-5000** \$140.00
- Sterile, Pk/1 (individually wrapped)
 - 5 mL Capacity (1 mL graduated)
 - Bx/400. **P200-5201S** \$36.95 ↓
 - Cs/10 Bx. **P200-521S** \$330.00
- Sterile, Pk/20
 - 5 mL Capacity (1 mL graduated)
 - Bx/500. **HS206371B** \$30.25
 - 5 mL Capacity (2 mL graduated)
 - Bx/500. **HS206372B** \$30.25




Pipettor Holders

- Bench Clamps, polypropylene, Heathrow
 - 1-Place, Clamp-Ette™
 - Blue. **HS2062A** \$12.90
 - Assorted (1 blue, 1 red & 1 purple). **HS2062B** \$27.30
 - 3-Place, Tri-Clamp-Ette™
 - Blue. **HS20162A** \$26.80
 - Red. **HS20162B** \$26.80
 - Purple. **HS20162C** \$26.80



- Magnetized
 - 2-Place, purple ABS plastic, MagPette™. **HS23503** \$16.20




- Stands
 - 4-Place, clear acrylic. **HS206204** \$35.25
 - 6-Place
 - Blue, opaque ABS plastic. **HS20613E** \$34.85
 - Clear, acrylic. **HS20620** \$45.65
 - Fluorescent Blue, transparent acrylic. **HS20613C** \$47.15
 - Green, opaque ABS plastic. **HS20613F** \$34.85

Continued



Wall-Mount
5-Place, TrippNT, USA. **50696** \$28.00↓



Reagent Reservoirs, v-bottom, hold 55 mL, Heathrow
Nonsterile, clear PVC, Pk/100. **HS20521A** \$36.30
Sterile, white polystyrene, Pk/50. **HS20521C** \$63.25
Individually sealed in a polyethylene bag and gamma radiated.



Polypropylene



Stainless Steel

Pitchers
Polypropylene, graduated (mL), tall form, tapered, Bel-Art
500 mL, 10 mL graduations, each. **289900000** \$5.40
1000 mL, 50 mL graduations, each. **289910000** \$7.70
2000 mL, 100 mL graduations, each. **289920000** \$8.80
3000 mL, 100 mL graduations, each. **289930000** \$12.20
Stainless Steel, graduated (mL/oz), tapered, Polar Ware
500 mL/16 oz, each. **T1062** \$53.75
1000 mL/32 oz, each. **T1063** \$69.00
2000 mL/64 oz, each. **T1064** \$86.50

Rulers, vinyl, 15 cm/6" long
Transparent. **4461** \$0.76
White. **4460** \$0.48

Rockers, Shakers & Mixers



BR1000



B3D1020



B3D1320

Rockers, Benchmark
Fixed Speed
2-D (seesaw), 12 x 12" platform w/flat mat
BenchBlotter™. **BR1000** \$379.00

3-D (nutating)
Small, 10.5 x 7.5" platform w/dimpled mat
MiniMixer™. **B3D1020** \$350.00
Large, 12 x 12" platform w/dimpled mat
BioMixer™. **B3D1320** \$379.00



BR2000



B3D2300

Variable Speed
• Speed is adjustable from 2-30 rpm
• Platform angle adjusts from 0-30°
• Cold room & incubator safe (4-65°C)
• Load capacity is 2 kg (4.5 lbs)
• 2-year warranty; maintenance-free brushless motors
2-D (seesaw), 14 x 12" platform w/flat mat
BenchRocker™ 2D. **BR2000** \$712.00
3-D (nutating), 12 x 12" platform w/dimpled and flat mats
BenchRocker™ 3D. **B3D2300** \$663.00



BT300



BT3000

Orbital Shakers, Benchmark
Small, 9.5 x 8" platform w/flat mat, 19 mm orbit
Orbi-Shaker™ JR. **BT300** \$800.00
Large, 14 x 12" platform w/flat mat, 19 mm orbit
Orbi-Shaker™. **BT3000** \$1,328.00



Vortex Mixer, 200-3200 rpm, BenchMixer™. **BV1000** \$218.00



Specimen Boxes (Pill Boxes), matchbox style, USA
2 1/8 x 1 1/8 x 5/8", Pk/72. **47** \$32.80
2 9/16 x 1 9/16 x 1", Pk/72. **49** \$41.10
3 1/8 x 2 1/16 x 1 1/4", Pk/72. **49XX** \$49.70

Continued