

## Illuminators

*Please call regarding availability of items not listed.*

### Fiber Optics, Dolan-Jenner

- Light Source, 150W EKE lamp, rheostat. **MI-150**.....\$343.00
- Gooseneck, Dual, chrome plated, 23" long. **EEG2823M**.....\$225.00
- Focusing Lens. **LH-759**.....\$38.00



### Fiber Optics, Schott

- Light Sources, 150W cold light
- Rheostat control, ACE
- EKE. **A20500** .....\$358.00
- EJA (*brighter*). **A20500.1** .....\$358.00
- Rheostat & Iris diaphragm, ACE I
- EKE. **A20520** .....\$471.00
- EJA (*brighter*). **A20520.1** .....\$471.00
- Goosenecks, semi-obedient
- Single, Chrome plated, 18" long. **A08400**.....\$141.00
- Dual, Chrome plated, 18" long. **A08500** .....\$209.00
- Focusing Lens. **A08080** .....\$41.40
- Polarizing Cap for A08080. **A08090** .....\$51.20

### Ringlights

- Universal, 47-62 mm ID, 33" long. **A08600** .....\$336.00
- Polarizer/Analyzer. **A08615** .....\$172.00
- Midi-Ringlight, 38-53 mm ID, 36" long. **A08630**.....\$301.00

### Filters

- Daylight. **A08074** .....\$41.35
- Blue. **A08070** .....\$41.35



MA263



MA264H

- LED Illuminator, 5W LED. MA264/LED** .....\$405.00

- Miniature Flood Light, 30W daylight bulb. MA263/100**.....\$96.00
- Freestanding illuminator has long-life bulb.

- Nicholas Illuminator, 6V/20W halogen. MA264H/100**.....\$289.00
- Very versatile unit. Variable transformer provides focused light.

## Magnifiers




Large diameter (4 1/8") glass lens is supported on a sturdy plastic, foldout stand for hands-free viewing. Base is graduated in both cm and inches for measuring objects.

- Folding magnifier (Linen Tester), 3X. ELP-800.03** .....\$10.10



A 4X magnifying lens is attached to fine-tipped, stainless steel forceps. Forceps double as a lens support, freeing viewer's hands. Folds for storage.

- Forceps magnifier. M15** .....\$4.00



Lenses are lightweight, water-clear, and shatterproof. Molded handles protect the lens and offer comfort and balance.

### Handle magnifiers, Ultra-Optix, USA

- 2" Lens, 3X. **SV-2P**.....\$4.15
- 3" Lens, 2.5X with 6X bifocal inset. **SV-3P**.....\$5.80
- 4" Lens, 2X with 6X bifocal inset. **SV-4P**.....\$9.80



Provides precise binocular viewing. Visor tilts out of the way when optical glass magnifiers are not needed. Adjustable headband is leather lined. Other magnifications and replacement lenses are available.

### Headband magnifier, Optivisor®, USA

- 1.75X at 14", DA3. **ELP-503.00**.....\$40.95
- 2.5X at 8", DA5. **ELP-505.00**.....\$40.95
- Light Attachment, Quasar™ LED. **ELP-558.00** .....\$24.95



View small details in low light with this compact, illuminated magnifier. Protective case keeps the lens scratch-free and serves as a handle when lens is extended.

- Lighted magnifier, 2 x 2" lens, 2X. SV-BPL**.....\$6.05



ELP-174.10



ELP-750.01

### Loupes

- Eye, plastic, 10X at 1", black. **ELP-174.10**.....\$2.60
- Triplet, 10X at 1", 18 mm, black. **ELP-750.01** .....\$16.00

Three optical glass lenses bonded together create a single achromatic (color-corrected) lens. Leather case provided.



- Stand magnifier, 2X glass lens w/2 alligator clips. HOL-163.00**....\$9.75
- Heavy base holds lens and 2 clamps.