

Illuminators

Please call regarding availability of items not listed.



- Fiber Optics, Schott-Fostec**
Light Sources, AC, 150W cold light
 Rheostat control, ACE
 DDL (*bright*). A20500.2.....\$314.00
 EKE (*brighter*). A20500.....\$314.00
 EJA (*brightest*). A20500.1.....\$314.00
 Rheostat & Iris diaphragm, ACE I
 DDL (*bright*). A20520.2.....\$412.00
 EKE (*brighter*). A20520.....\$412.00
 EJA (*brightest*). A20520.1.....\$412.00
Light Sources, DC Regulated, 150W cold light
 Rheostat control, DCR III
 DDL (*bright*). A20800.2.....\$576.00
 EKE (*brighter*). A20800.....\$576.00
 EJA (*brightest*). A20800.1.....\$576.00
 Rheostat & Iris diaphragm, DCR III I
 DDL (*bright*). A20810.2.....\$684.00
 EKE (*brighter*). A20810.....\$684.00
 EJA (*brightest*). A20810.1.....\$684.00
Goosenecks, semi-obedient
 Single, Chrome plated, 18" long. A08400.....\$124.00
 Single, Black tubing, 30" long. A08410.....\$126.00
 Dual, Chrome plated, 18" long. A08500.....\$182.00
 Focusing Lens. A08080.....\$36.00
 Polarizing Cap for A08080. A08090.....\$38.00
Ringlights
 Universal, 47-62 mm ID, 33" long. A08600.....\$294.00
 Polarizer/Analyzer. A08615.....\$151.00
 Midi-Ringlight, 38-53 mm ID, 36" long. A08630.....\$264.00
 Polarizer/Analyzer. A08632.....\$152.00
Filters
 Daylight. A08074.....\$34.00
 Blue. A08070.....\$32.00
 Green. A08071.....\$32.00
 Yellow. A08072.....\$32.00



- Miniature Flood Light, 30W daylight bulb. MA263/100.....\$96.00**
 Freestanding illuminator has long-life bulb.
- Nicholas Illuminator, 6V/20W halogen. MA264H/100.....\$289.00**
 Very versatile unit. Variable transformer provides focused light.

Magnifiers



Large diameter glass lens is supported on a sturdy plastic, foldout stand for hands-free viewing. Base is graduated in both cm and inches for measuring objects.

- Folding Magnifier (Linen Tester), 4 1/8" lens, 3X.....\$10.10**



A 4X magnifying lens is attached to fine-tipped, stainless steel forceps. Forceps double as a lens support, freeing viewer's hands. Folds for storage.

- Forceps Magnifier. M15.....\$4.00**



Lenses are lightweight, water-clear, and shatterproof. Molded handles protect the lens and offer comfort and balance.

- Handle Magnifiers, Ultra-Optix**
 2" Lens, 3X. SV-2P.....\$4.15
 3" Lens, 2.5X with 6X bifocal inset. SV-3P.....\$5.80
 4" Lens, 2X with 6X bifocal inset. SV-4P.....\$9.80
 2 x 4" Lens, 2X with 6X bifocal inset. SV-XP.....\$9.50



Provides precise binocular viewing. Visor tilts out of the way when optical glass magnifiers are not needed. Adjustable headband is leather lined. Other magnifications and replacement lenses are available.

- Headband Magnifier, Optivisor®**
 1.75X at 14", DA3. ELP-503.00.....\$39.95
 2.5X at 8", DA5. ELP-505.00.....\$39.95



View small details in low light with this compact, illuminated magnifier. Protective case keeps the lens scratch-free and serves as a handle when lens is extended.

- Lighted Magnifier, 2 x 2" lens, 2X. SV-BPL.....\$6.05**



- Loupes**
 Eye, plastic, 10X at 1", black. ELP-174.10.....\$2.60
 Triplet, 10X at 1", 18 mm, black. ELP-750.01.....\$16.00
 Three optical glass lenses bonded together create a single achromatic (color-corrected) lens. Leather case provided.



- Stand Magnifier, 2X glass lens with 2 alligator clips.....\$9.75**
 Heavy base holds lens and 2 clamps.