

Tubes



- Centrifuge, conical, polypropylene, sterile
15 mL, sterile bag/50, Cs/500. **HS4426B** \$83.20
- 50 mL, max RCF is 9,500 x g
Natural, sterile bag/50, Cs/500. **HS4427B** \$114.00
- Black, sterile bag/50, Cs/500. **HS4429** \$165.00



- Centrifuge w/Attached Cap, polypropylene, max RCF is 5,000 x g
5 mL, 0.5 mL graduations, green, Pk/250. **T2079** \$75.95



Polystyrene with dual-position caps.

- Culture, 12 x 75 mm (5 mL), sterile, Simport
Cs/500, sterile bag/25. **T415-2** \$63.00
- Cs/1000, sterile bag/125. **T415-3** \$94.00




Natural



Black

- Microcentrifuge, Economy, polypropylene (PP), natural color
0.5 mL, 0.1 mL graduations, natural, Pk/1000. **HS4422** \$32.20
- 1.5 mL, 0.25 mL graduations
Natural, max RCF is 17,000 x g
Pk/500. **HS4323** \$15.05
- Pk/1000. **HS4325** \$29.95
- Black, max RCF is 18,000 x g, Pk/500. **HS4323K** \$17.40



0.6 mL

1.5 mL

2.0 mL

Tubes will withstand centrifugal forces over 30,000 x g and temperatures from -80° to 121°C. Flat, textured caps are pierceable.

- Microcentrifuge, Premium, PP, G-Tube™, Pk/500, USA
Amber color protects light-sensitive specimens.
- 0.6 mL \$24.65
- Amber. **4048** Green. **4044** Natural. **4040**
- Blue. **4043** Lavender. **4047** Yellow. **4042**
- Mixed (blue, green, orange, pink, & yellow). **4049** \$29.70
- 1.5 mL, 0.5 mL graduations, frosted area on side \$24.65
- Amber. **4038** Green. **4034** Natural. **4030**
- Blue. **4033** Lavender. **4037** Yellow. **4032**
- Mixed (blue, green, orange, pink, & yellow). **4039** \$29.70

Continued

- 2.0 mL, 0.5 mL graduations, frosted area on side \$24.10
- Amber. **4058** Green. **4054** Natural. **4050**
- Blue. **4053** Lavender. **4057** Yellow. **4052**
- Mixed (blue, green, orange, pink, & yellow). **4059** \$31.55

- Microcentrifuge, Siliconized, G-Tube™, Pk/500, USA
1.5 mL, natural color. **4165SL** \$46.00
(Siliconized products should NOT be autoclaved)




Securely close lids on 0.5 and 1.5 mL microcentrifuge tubes with these autoclavable, polypropylene holders. Extract tubes from centrifuges, dry baths, racks, etc. using large tabs. Assorted colors.

- Microcentrifuge Tube Holders, Pk/100. **HS4350** \$21.35

- PCR, 0.2 mL thin-wall, polypropylene, USA
Strips (*caps are sold separately; available in natural only*)
- 8-Tube Strips, Pk/125 \$82.25
- Blue. **5010-4** Natural. **5010-1**
- 8-Tube Cap Strips, domed, Pk/300. **5030-1** \$80.55
- 12-Tube Strips, Pk/100 \$100.00
- Blue. **5020-4** Natural. **5020-1**
- 12-Tube Cap Strips, domed, Pk/200. **5040-1** \$80.55

Tube Racks



- Floating, PVC Foam, Bel-Art, USA
- Circular. *Recommended for ice buckets*
- 44-Hole, 1.5 mL tubes, 22 cm dia. x 1 cm \$28.05
- Black/**188480012** Green/**188480014** Red/**188480013**
- 30-Hole, 14/15 mL tubes, 22 cm dia. x 1.7 cm \$29.95
- Black/**188480022** Green/**188480024** Red/**188480023**
- Rectangular. *Recommended for ice pans*
- 36-Hole, 1.5 mL tubes, 23.5 x 11 x 1 cm \$28.05
- Black/**188480112** Green/**188480114** Red/**188480113**
- 24-Hole, 14/15 mL tubes, 23.5 x 11 x 1.7 cm \$29.95
- Black/**188480122** Green/**188480124** Red/**188480123**



- Floating, HDPE Foam w/Handle, Pk/5
- 18-Hole, 0.2/0.5/1.5-2.0 mL tubes, blue circle. **HS2166** \$13.60
- 24-Hole, 1.5-2.0 mL tubes, yellow rectangle. **HS2165C** \$13.40
- 8-Hole, 15 mL tubes, blue parallelogram. **HS2165B** \$11.80
- 4-Hole, 50 mL tubes, green diamond. **HS2165A** \$13.20

Continued