

HS29025



HS2345

Microtube Racks

Polypropylene

- 80-Well (16 x 5), molded-in coordinates, 22.5 x 6.7 x 2.8 cm
Holds 1.5-2.0 mL tubes in conical wells.
- Blue, Pk/5. **HS29025B** \$23.75
- Pink, Pk/5. **HS29025D** \$23.75
- 96-Well reversible (12 x 8), snap-on cover, 24.6 x 12.1 x 5 cm
Holds 1.5-2.0 mL tubes on one side or 0.5 mL on reverse.
- Blue, each. **HS2345C** \$13.50
- Green, each. **HS2345D** \$13.50
- Pink, each. **HS2345E** \$13.50

Holds 1.5-2.0 mL tubes
in a two-piece rack.



- Polystyrene Foam (styrofoam)
- 50-Well, Pk/10. **HS4283P** \$36.90
- 100-Well, Pk/10. **HS4282** \$46.45



- PCR Reversible racks, assorted colors, Pk/5. **HS2344A** \$47.95
- Holds 0.2, 0.5 and 1.5 mL tubes.



- Peg racks, polypropylene
- 96-Place, 10-13 mm tubes, blue, each. **HS24312B** \$11.80
- 50-Place, 14-17 mm tubes, blue, each. **HS24311B** \$10.70



- Refrigerator Racks, two-piece polystyrene
- 50-Well, 1.5 mL tubes, each. **HS28230A** \$35.35
- 15-Well, 15/50 mL tubes, each. **HS28230B** \$35.35



- Slant racks, 40-place, polypropylene, Bel-Art
- Holds 13-16 mm Tubes, each. **188900016** \$46.75
- Holds 16-20 mm Tubes, each. **188900020** \$46.75




Tilt your test tube rack 25°, 35°, or 45°. Create agar slant tubes, or adjust the angle of tube racks when pipetting to reduce hand strain.

- Tube Rack Slanter. **188900100** \$15.95
- Platform size: 19.6 cm (W) x 12.7 cm (D)



- Stackable, 72-place, polypropylene, Bel-Art
- Holds 1.5 mL Tubes, each. **188410000** \$32.80



- SmoothRack™, 72-place, polypropylene, blue, USA
- Holds 10-17 mm Tubes, each. **SR-B** \$14.75



Remains stable
in water bath.

- Snap-Together racks, polypropylene, blue
- 90-Place, 13 mm tubes, each. **HS24305A** \$12.70
- 60-Place, 16 mm tubes, each. **HS24305B** \$12.70
- 40-Place, 19 mm tubes, each. **HS24305C** \$12.70