

MAINTENANCE & REPLACEMENT PARTS

Bulbs/Lamps



Please call regarding availability of bulbs not listed.

ANSI codes and manufacturer's part numbers are cross-referenced (in italics and parenthesis) when available.

Tungsten Halogen & Incandescent

* Recommend Bulb Removal/Installation Tool (shown after listing)

2.5V/0.35A Bullet. 252	\$4.00
<i>(634, 8000155)</i>	
6V/0.2A Bullet (for teaching head). GE328	\$2.00
<i>(79098, 31-31-50, 8000168)</i>	
6V/0.2A Bullet w/Lens (for multi-TH), Nikon. MXA20101	\$11.00
<i>(79186)</i>	
6V/1.2A Horizontal Filament, Hosobuchi. OP2101	\$9.00
<i>(MA560, 8000369)</i>	
6V/1.2A Vertical Filament, Hosobuchi. OP2101/U	\$9.00
<i>(MA561, 8000370)</i>	
6V/2A SC Bayonet, Frosted, Ushio. 8-C201	\$7.10
<i>(8000301)</i>	
6V/2A SC Bayonet, Round. OP2106	\$6.00
<i>(8-M102, 8000305)</i>	
6V/3A (6V/18W) Frosted, ½ Silver, Swift. MA737	\$10.00
<i>(8000350)</i>	
6V/4.5A DC Prefocus*, Ushio. 1468	\$5.00
<i>(581, 8000067)</i>	
6.4V/2.63A DC Bayonet. S8.1130	\$2.30
6.5V/2.75A (6.5V/15W) DC Bayonet. 1493	\$2.75
<i>(360, 700-034, 8-B102, MA402, 21138, 8000070)</i>	
6.5V/2.75A DC Prefocus*. GE1460	\$9.60
<i>(MA760, 1033, 8000065)</i>	
6.5V/2.75A DC Prefocus*, Ushio. 1460X	\$5.50
<i>(31-31-75, 31-31-76, 31-01-05, 8000066)</i>	
6.5V/2.75A DC Prefocus*. GE1630	\$9.30
<i>(8000092)</i>	
6.5V/2.75A DC Prefocus*. 1631X	\$8.00
<i>(31-31-85, 31-31-86, 31-31-79, 8000093)</i>	
6.5V/2.75A DC Prefocus*. GE1631X	\$13.50
<i>(31-31-85, 31-31-86, 31-31-79, 8000093)</i>	
6.5V/2.75A DC Prefocus*, Ushio. 1649	\$9.90
<i>(31-31-42, 31-31-43, 8000100)</i>	
6V/5W Screw, Ushio. 8000302	\$10.00
<i>(8-G101)</i>	
6V/9W MR11, Philips. 13529	\$21.50
6V/10W (6V/1.67A) Halogen, G4. GE787	\$7.50
<i>(31-31-77-01, 31-31-77-10, MA570, 76342, 8000254)</i>	
6V/10W Halogen, G4, Osram. 64223	\$17.00
<i>(38-61-08, 91-01-54, 5972, 1000857)</i>	

6V/10W Halogen, G4, Osram. 64225	\$4.50
<i>(ESA/FHD, 76342, 050-550, 050-556, 700-065, 8-C403, 907-090, 624, 655C, 7387, 3369, 1000531)</i>	
6V/10W Halogen, G4, Philips. 5972	\$7.50
<i>(38-61-08, 91-01-54, 64223, 1000857)</i>	
6V/10W Halogen, G4, Philips. 7387	\$5.20
<i>(ESA/FHD, 76342, 050-550, 050-556, 700-065, 8-C403, 907-090, 624, 655C, 7387, 3369, 1000531)</i>	
6V/10W SC Bayonet, round. 5001	\$1.83
<i>(3040-0260, 6814P)</i>	
6V/10W Screw, Hosobuchi. OP2103	\$8.00
<i>(8-G102, 8000303)</i>	
6V/15W MR11, Philips. 13528	\$22.00
<i>(8-C408, 630)</i>	
6V/15W Bulb and Socket, Ushio. 8-C103	\$18.00
<i>(8000300)</i>	
6V/15W DC Bayonet, Flanged, Oblong. 77903	\$10.00
<i>(OP2130, 8000240)</i>	
6V/15W DC Bayonet, Flanged, Nikon. 77910	\$13.50
<i>(8000242)</i>	
6V/15W DC Bayonet, Oblong, Ushio. 8018	\$7.40
<i>(3800-29-9020, 76906, 166-324, 500-012, 8000269)</i>	
6V/15W DC Flanged. 76Z	\$17.90
<i>(3800-18-1730, 38-01-76, 8000170)</i>	
6V/15W DC Flanged. 177Z	\$9.25
<i>(3800-18-1740, 38-01-77, OP2131, 8000172)</i>	
6V/15W (6V/3A) DC Bayonet, Round. OP2114	\$6.25
<i>(76907, 8000233)</i>	
6V/20W Halogen, G4. GE778	\$8.90
<i>(31-31-81-01, 31-31-81-10, 569, 31-31-21, 8000239)</i>	
6V/20W (6V/3.33A) Halogen, G4. GE788	\$7.90
<i>(MA260/05, 1121)</i>	
6V/20W Halogen, G4, Osram. 64250 HLX	\$4.40
<i>(ESB/FHE, 1120, 362-658, 3801-20-1740D, 11143, 500-245, 89500, 79099, 8-C405, MA780, 7388, M420, MA07, 687.95, 3256, 12-561-BLB1, MA311, 1000532)</i>	
6V/20W Halogen with Lead, Narva. 55133	\$33.00
6V/20W MR11, Ushio. JCR/M6V20W	\$18.40
<i>(83314)</i>	
6V/25W Flanged, G4. 39-01-53	\$19.00
<i>(8000174)</i>	
6V/25W Halogen with Lead, Narva. 55147	\$33.00
<i>(33-01-53)</i>	
6V/30W (6V/5A) Bayonet, Oblong, Hosobuchi. OP2118	\$6.70
<i>(MA738, 77904, 8-C101, 8000241)</i>	
6V/30W Bulb and Socket, Ushio. 8-C102	\$17.50
<i>(8000299)</i>	
6V/30W (6V/5A) Candelabra, Osram. 8100	\$18.75
<i>(500-011, 8000277)</i>	
6V/30W DC Flanged, Ushio. 8000175	\$20.00
<i>(39-01-58)</i>	
6V/30W Halogen, G4, Philips. 5761	\$7.40
<i>(84246, 8-C410, MA326, 3256-30W, LB-30WH, SE06.0215, 1000867)</i>	
6V/30W Halogen, G4, Ushio. 1000867	\$11.90
<i>(84246, 8-C410, MA326, 3256-30W, LB-30WH, SE06.0215, 5761)</i>	
6V/30W Halogen with Lead, Narva. 55168	\$39.00
6V/35W Halogen, G4, vertical filament, Osram. 64275	\$7.50

Continued