

Illuminators

Please call regarding availability of items not listed.

Fiber Optics, Dolan-Jenner

- Light Source, 150W EKE lamp, rheostat. **MI-150**.....\$343.00
- Gooseneck, Dual, chrome plated, 23" long. **EEG2823M**.....\$225.00
- Focusing Lens. **LH-759**.....\$38.00



Fiber Optics, Schott

- Light Sources, 150W cold light
- Rheostat control, ACE
- EKE. **A20500**\$358.00
- EJA (*brighter*). **A20500.1**\$358.00
- Rheostat & Iris diaphragm, ACE I
- EKE. **A20520**\$471.00
- EJA (*brighter*). **A20520.1**\$471.00
- Goosenecks, semi-obedient
- Single, Chrome plated, 18" long. **A08400**.....\$141.00
- Dual, Chrome plated, 18" long. **A08500**\$209.00
- Focusing Lens. **A08080**\$41.40
- Polarizing Cap for A08080. **A08090**\$51.20

Ringlights

- Universal, 47-62 mm ID, 33" long. **A08600**\$336.00
- Polarizer/Analyzer. **A08615**\$172.00
- Midi-Ringlight, 38-53 mm ID, 36" long. **A08630**.....\$301.00

Filters

- Daylight. **A08074**\$41.35
- Blue. **A08070**\$41.35



MA263



MA264H

- LED Illuminator, 5W LED. MA264/LED**\$405.00

- Miniature Flood Light, 30W daylight bulb. MA263/100**.....\$96.00
- Freestanding illuminator has long-life bulb.

- Nicholas Illuminator, 6V/20W halogen. MA264H/100**.....\$289.00
- Very versatile unit. Variable transformer provides focused light.

Magnifiers




Large diameter (4 1/8") glass lens is supported on a sturdy plastic, foldout stand for hands-free viewing. Base is graduated in both cm and inches for measuring objects.

- Folding magnifier (Linen Tester), 3X. ELP-800.03**\$10.10



A 4X magnifying lens is attached to fine-tipped, stainless steel forceps. Forceps double as a lens support, freeing viewer's hands. Folds for storage.

- Forces magnifier. M15**\$4.00



Lenses are lightweight, water-clear, and shatterproof. Molded handles protect the lens and offer comfort and balance.

Handle magnifiers, Ultra-Optix, USA

- 2" Lens, 3X. **SV-2P**.....\$4.15
- 3" Lens, 2.5X with 6X bifocal inset. **SV-3P**.....\$5.80
- 4" Lens, 2X with 6X bifocal inset. **SV-4P**.....\$9.80



Provides precise binocular viewing. Visor tilts out of the way when optical glass magnifiers are not needed. Adjustable headband is leather lined. Other magnifications and replacement lenses are available.

Headband magnifier, Optivisor®, USA

- 1.75X at 14", DA3. **ELP-503.00**.....\$40.95
- 2.5X at 8", DA5. **ELP-505.00**.....\$40.95
- Light Attachment, Quasar™ LED. **ELP-558.00**\$24.95



View small details in low light with this compact, illuminated magnifier. Protective case keeps the lens scratch-free and serves as a handle when lens is extended.

- Lighted magnifier, 2 x 2" lens, 2X. SV-BPL**.....\$6.05



ELP-174.10



ELP-750.01

Loupes

- Eye, plastic, 10X at 1", black. **ELP-174.10**.....\$2.60
- Triplet, 10X at 1", 18 mm, black. **ELP-750.01**\$16.00

Three optical glass lenses bonded together create a single achromatic (color-corrected) lens. Leather case provided.



- Stand magnifier, 2X glass lens w/2 alligator clips. HOL-163.00**....\$9.75
- Heavy base holds lens and 2 clamps.

Filters, glass



29.8 mm Diameter	
Blue, Clear. MA289/05	\$30.00
Blue, Frosted. MA289	\$30.00
Diffusion (frosted). MA288	\$30.00
Green, clear. MA290/05	\$30.00
32 mm Diameter	
Blue, Clear. 923-160	\$10.00
Blue, Frosted. 3236	\$7.50
Diffusion (frosted). 3235-W	\$7.50
Green, clear. CM-0321-G	\$7.00
40 mm Diameter	
Blue, frosted. MA563/05	\$45.00
Diffusion (frosted). MA566	\$45.00
45 mm Diameter	
Blue, clear. 14671	\$10.40
Color Balancing, NCB-11, Nikon. MBN11710	\$77.00
Diffusion (frosted). 14673	\$35.00
Didymium (H&E), BG-20, Olympus. 9-U116	\$52.00
Green Interference, 546/60, Nikon. MBN11200	\$52.00
Neutral Density, Nikon	
50% Transmittance, ND2A. MBN21802	\$42.00
12.5% Transmittance, ND8A. MBN21808	\$42.00

Fluorescence



Filter Holders, empty, USA	
For Leica Ploempaks. 91004	\$350.00
For Nikon Microscopes with 18 mm Filter Ports. 91000	\$150.00
For Nikon Microscopes with 25 mm Filter Ports. 91001	\$350.00
For Olympus BH2 Microscopes. 91002	\$275.00
For Olympus BX Microscopes. 91018	\$385.00
Filter Sets, Chroma, USA.....	Access www.chroma.com
	and call us to place an order



Confirm patient identification and preview the slide before examining microscopically. This macro inspection magnifier mounts onto the microscope's arm, where the lens angle can be adjusted using the double-ball joint arm. Large, 3 x 1¼ " lens permits viewing of entire 3 x 1" glass slide.

Magnifier, "Martin's Macro Magnifier", USA. 15401	\$12.50
--	---------



Nosepiece Dust Plug, RMS thread (20.32 mm Dia.). 15420	\$0.50
---	--------

Object Markers, RMS thread (20.32 mm Dia.)

Mark slides with a 1.8 mm ring for later referral.

Black. 14361	\$172.00
Purple. 14381	\$172.00
Red. 14371	\$172.00
Refill Ink, 10 mL, Shachihata	
Black. 22112	\$7.00
Purple. 22115	\$7.00
Red. 22111	\$7.00

Adapters

RMS/45 to CFI60, 15 mm long. 14111	\$43.00
RMS ID to 25 mm OD Thread, brass. 14115	\$52.00

Microscopes with Nikon CFI60 infinity-corrected (∞) objectives require a 14111 adapter.
Microscopes with Leica 25 mm thread/45 mm parfocal length objectives require a 14115 adapter.
Microscopes with 160 mm tube length objectives do not require an adapter.



Objectives, Jenco, Meiji, National, Nikon, Olympus

160 mm tube length, DIN

Achromats, Motic/National

4X, N.A. 0.10. 704-160	\$30.00
10X, N.A. 0.25. 710-160	\$31.00
40X, N.A. 0.65. 740-160	\$40.00
100X Oil, N.A. 1.25. 799-160	\$65.00

Super High Contrast Achromats, National

4X, N.A. 0.10. 704-160ASC	\$26.00
10X, N.A. 0.25. 710-160ASC	\$34.00
40X, N.A. 0.65. 740-160ASC	\$62.00
100X Oil, N.A. 1.25. 799-160ASC	\$95.00

Semi-Plan Achromats, National (unless noted)

2X, N.A. 0.05, Accu-Scope. 3193SP-00	\$70.00
4X, N.A. 0.10. 704-160SP	\$25.00
10X, N.A. 0.25. 710-160SP	\$48.00
40X, N.A. 0.65. 740-160SP	\$95.00
50X Oil, N.A. 0.95, Meiji. MA694	\$280.00
60X, N.A. 0.85. MA645	\$300.00
100X Oil, N.A. 1.25. 799-160SP	\$150.00

Plan Achromats, National (unless noted)

2.5X, N.A. 0.08, Meiji. MA432	\$270.00
4X, N.A. 0.10. 704-160P	\$31.00
10X, N.A. 0.25. 710-160P	\$52.00
20X, N.A. 0.40. 720-160P	\$115.00
40X, N.A. 0.65. 740-160P	\$99.00
50X Oil, N.A. 0.95, Accu-Scope. 3190PL-00	\$315.00
60X, N.A. 0.85. 760-160P	\$185.00
100X Oil, N.A. 1.25. 799-160P	\$160.00