

- 120V/300W MR16, stippled, oval 2-pin, Osram. **ELH** ..... \$11.75  
(64660, 1000321)
- 125V/300W MR16, Long Life, faceted, Eiko. **ELH-5** ..... \$17.00  
(1003020)
- 130V/50W Halogen. **RN13050** ..... \$40.00

**Metal Halide**

- 60V/24W MR11, 300 hrs, Welch Allyn. **09500-U**..... \$270.00  
(AL-1824, 8-61257)
- 120W Reflector, 1500 hrs, Lumen Dynamics. **P012-63000** .. \$575.00  
(OEX-BULB)

**Mercury Short Arc**

- 50W Mercury Arc, 100 Hrs, Osram. **HBO50W/AC/L1**..... \$115.00  
(78564, 78588, 8-B199, 38-16-19, 414140, 2055,  
3040-0090, 500-137, SHL-50AC1, BA-001, 5000487)
- 50W Mercury Arc, 100 Hrs, Osram. **HBO50W/AC/L2**..... \$115.00  
(31-31-65, SHL-50AC2, BA-002, 5000857)
- 50W Mercury Arc, 200 Hrs, Osram. **HBO50W/3** ..... \$140.00  
(BA-003)
- 100W Mercury Arc, 200 Hrs, Ushio. **USH-102D** ..... \$109.00  
(Burn in vertical position only)  
(78589, 8-B192U, 500-138, 3800-18-4060/D, 2087-001,  
3040-0380, BA-005A, HBO100W/2, M/5000/133, 5000274)
- 100W Mercury Arc, 200 Hrs, Ushio. **USH-102DH** ..... \$109.00  
(Burn in vertical or horizontal position)  
(78589, 8-B192U, 500-138, 3800-18-4060/D, 2087-001,  
3040-0380, BA-005A, HBO100W/2, M/5000/133, 5000274)
- 100W Mercury Arc, 300 Hrs, Ushio. **USH-103D**..... \$123.00  
(78590, 8-B1920, 380301-9350, USH-103OL, 5001329)
- 100W Mercury Arc, 200 Hrs, Osram. **HBO100W/2** ..... \$109.00  
(78589, 8-B192U, 500-138, 3800-18-4060/D, M/5000/133  
2087-001, 3040-0380, BA-005, 3254, USH-102DH)
- 100W Mercury Arc, 300 Hrs, Osram. **HBO103W/2** ..... \$123.00  
(Not compatible with all Hg power supplies)  
(78590, 8-B1920, 380301-9350, 500-321)
- 200W Mercury Arc, 1000 hrs, Ushio. **USH-200DP**..... \$135.80  
(HBO200W/2, 31-01-01, 649, 860014, 78580, 8-B190,  
3800-18-4070, 38-16-15, 050-320, MA885,  
91-00-50, 505-620, 78590, 5000276)

**Xenon Short Arc**

- 75W Xenon Arc, 400 hrs, Ushio. **UXL-75XE** ..... \$195.00  
(XBO75W/2, 8-B198, 910881-2013, 38-16-02, 500-139,  
91-01-42, 3800-79-9190, 78682, 620-125, 903-050,  
BA-013, 5000371)
- 75W Xenon Arc, 400 hrs, Osram. **XBO75W/2** ..... \$155.00  
(UXL-75XE, 8-B198, 910881-2013, 38-16-02, 500-139,  
91-01-42, 3800-79-9190, 78682, 620-125, 903-050,  
BA-013, 5000371)
- 100W Xenon Arc, 500 hrs, Osram. **XBO100WOFR**..... \$300.00  
(BA-014A)
- 150W Xenon Arc, Short Gap. **XBO150W/CROFR** ..... \$260.00  
(8-B195, BA-018, UXL-150S-0, 3801-22-2450,  
500-008, 5000326)

- 150W Xenon Arc, Standard Gap. **XBO150W/1OFR** ..... \$214.00  
(8-B196, BA-016)
- 150W Xenon Arc, Long Life. **UXL-S150M-O**..... \$355.00
- 450W Xenon Arc, 2000 hrs, Osram. **XBO450W/2OFR**..... \$460.00  
(UXL-450S-O, 5000355)
- 450W Xenon Arc, 2000 hrs, Ushio. **UXL-450S-O** ..... \$535.00  
(XBO450W/2OFR, 5000355)

**Fluorescent**

- 4W Fluorescent, Tubular, Eiko. **F4T5/D** ..... \$7.00  
(31-31-03, MA309/05, 3000108)
- 5W Fluorescent, G23, Sylvania. **CF5DS/841** ..... \$5.00  
(800-138, MA68, 48920-81, 20303)
- 9W Fluorescent, G23, Sylvania. **CF9DS/841** ..... \$4.95  
(20305)
- 22W Fluorescent, 8" OD Circle, Sylvania. **FC8T9/D/RS** ..... \$19.00  
(29980)
- 32W Fluorescent, 12" OD Circle, Sylvania. **FC12T9/D/RS** ..... \$15.90  
(20030)

**Germicidal**

- 15W Germicidal, 2-pin, Sylvania. **G15T8** ..... \$17.00  
(21612, 3000007)
- 30W Germicidal (UV-C), 2-pin, Sylvania. **G30T8** ..... \$20.00  
(23112, 00204, 11080, 3000009)
- 39W Germicidal, single-pin, Sylvania. **G36T5L**..... \$28.80  
(23443, 00206)
- 65W Germicidal, single-pin, 155.4 cm, Higuchi. **G64T5L** ..... \$83.00

**Heat (Infrared)**

- 250W, R40, shatter resistant, GE. **250BR40/HI/TC** ..... \$14.00

\* Recommend Bulb Removal/Installation Tool (shown below)



Remove and install prefocused/flanged bulbs safely with this simple tool. Push rubber-lined tube directly onto top of bulb and turn clockwise to install, or counterclockwise to remove.

- Bulb Removal/Installation Tool. 3283** ..... \$16.00

- Compressed Air, 10 oz Can. 77001** ..... \$11.00  
Ultra pure duster is nonflammable.



Continued