



- Dust Covers, Vinyl, translucent, sewn-on binding seam, Vikem®
 Bag Style (bottom x top x height)
 16 x 11 x 18". 243030000\$11.30
 18 x 15 x 22". 243040000\$13.90
- England Field Finder, 1 mm thick. KR-834\$159.00

- Eyepiece Pointer, wire. 14301\$2.00
- 

- 
- Eyetebe Sleeve, 23.2 mm increased to 30 mm. 14320\$35.00

- Filters, glass
- 
- 29.8 mm Diameter
 Blue, Clear. MA289/05\$30.00
 Blue, Frosted. MA289\$30.00
 Diffusion (frosted). MA288\$30.00
 Green, clear. MA290/05\$30.00
- 32 mm Diameter
 Blue, Clear. 923-160\$8.00
 Blue, Frosted. 3236\$7.50
 Diffusion (frosted). 3235-W\$7.50
 Green, clear. CM-0321-G\$7.00
- 40 mm Diameter
 Blue, frosted. MA563/05\$45.00
 Diffusion (frosted). MA566\$45.00
- 45 mm Diameter
 Blue, clear. SG060727\$10.40
 Color Balancing, NCB-11, Nikon. MBN11710\$74.00
 Diffusion (frosted). 14673\$35.00
 Didymium (H&E), BG-20, Olympus. 9-B302\$51.00
 Green Interference, 546/60, Nikon. MBN11200\$51.00
 Neutral Density, Nikon
 12.5% Transmittance, ND8A. MBN21808\$37.00
 6.25% Transmittance, ND16A. MBN21816\$37.00

- Fluorescence
- 
- Filter Holders, empty
 For Leica Ploempaks. 91004\$350.00
 For Nikon Microscopes with 18 mm Filter Ports. 91000\$150.00
 For Nikon Microscopes with 25 mm Filter Ports. 91001\$350.00
 For Olympus BH2 Microscopes. 91002\$275.00
 For Olympus BX Microscopes. 91018\$385.00
 Filter Sets, ChromaAccess www.chroma.com
 and call us to place an order




Hoffman Modulation Contrast
 Modulation Optics..... Access www.modulationoptics.com
 and call us to place an order

Produce vivid, 3-dimensional images when using plastic or glass chambers. Systems are available for Leica, Nikon, Olympus, or Zeiss microscopes.



Magnifier, "Martin's Macro Magnifier". 15401\$12.50



Nosepiece Dust Plug, RMS thread (20.32 mm Dia.). 15420\$0.50

- Object Markers, RMS thread
 Mark slides with a 1.8 mm ring for later referral.
 Black. 14361\$166.00
 Purple. 14381\$166.00
 Red. 14371\$166.00
 Refill Ink, 10 mL, Shachihata
 Black. 22112\$7.00
 Purple. 22115\$7.00
 Red. 22111\$7.00
 Adapters
 RMS/45 to CFI60, 15 mm long. 14111\$43.00
 RMS ID to 25 mm OD Thread, brass. 14115\$52.00

*Microscopes with Nikon CFI60 infinity-corrected (∞) objectives require a 14111 adapter.
 Microscopes with Leica 25 mm thread/45 mm parfocal length objectives require a 14115 adapter.
 Microscopes with 160 mm tube length objectives do not require an adapter.*