

**Filters, glass**



29.8 mm Diameter	
Blue, Clear. <b>MA289/05</b> .....	\$30.00
Blue, Frosted. <b>MA289</b> .....	\$30.00
Diffusion (frosted). <b>MA288</b> .....	\$30.00
Green, clear. <b>MA290/05</b> .....	\$30.00
32 mm Diameter	
Blue, Clear. <b>923-160</b> .....	\$10.00
Blue, Frosted. <b>3236</b> .....	\$7.50
Diffusion (frosted). <b>3235-W</b> .....	\$7.50
Green, clear. <b>CM-0321-G</b> .....	\$7.00
40 mm Diameter	
Blue, frosted. <b>MA563/05</b> .....	\$45.00
Diffusion (frosted). <b>MA566</b> .....	\$45.00
45 mm Diameter	
Blue, clear. <b>14671</b> .....	\$10.40
Color Balancing, NCB-11, Nikon. <b>MBN11710</b> .....	\$77.00
Diffusion (frosted). <b>14673</b> .....	\$35.00
Didymium (H&E), BG-20, Olympus. <b>9-U116</b> .....	\$52.00
Green Interference, 546/60, Nikon. <b>MBN11200</b> .....	\$52.00
Neutral Density, Nikon	
50% Transmittance, ND2A. <b>MBN21802</b> .....	\$42.00
12.5% Transmittance, ND8A. <b>MBN21808</b> .....	\$42.00

**Fluorescence**



Filter Holders, empty, USA	
For Leica Ploempaks. <b>91004</b> .....	\$350.00
For Nikon Microscopes with 18 mm Filter Ports. <b>91000</b> .....	\$150.00
For Nikon Microscopes with 25 mm Filter Ports. <b>91001</b> .....	\$350.00
For Olympus BH2 Microscopes. <b>91002</b> .....	\$275.00
For Olympus BX Microscopes. <b>91018</b> .....	\$385.00
Filter Sets, Chroma, USA.....	Access <a href="http://www.chroma.com">www.chroma.com</a>
	and call us to place an order



Confirm patient identification and preview the slide before examining microscopically. This macro inspection magnifier mounts onto the microscope's arm, where the lens angle can be adjusted using the double-ball joint arm. Large, 3 x 1¼ " lens permits viewing of entire 3 x 1" glass slide.

<b>Magnifier, "Martin's Macro Magnifier", USA. 15401</b> .....	\$12.50
--	---------



<b>Nosepiece Dust Plug, RMS thread (20.32 mm Dia.). 15420</b> .....	\$0.50
---	--------

**Object Markers, RMS thread (20.32 mm Dia.)**


Mark slides with a 1.8 mm ring for later referral.

Black. <b>14361</b> .....	\$172.00
Purple. <b>14381</b> .....	\$172.00
Red. <b>14371</b> .....	\$172.00
Refill Ink, 10 mL, Shachihata	
Black. <b>22112</b> .....	\$7.00
Purple. <b>22115</b> .....	\$7.00
Red. <b>22111</b> .....	\$7.00

**Adapters**

RMS/45 to CFI60, 15 mm long. <b>14111</b> .....	\$43.00
RMS ID to 25 mm OD Thread, brass. <b>14115</b> .....	\$52.00

*Microscopes with Nikon CFI60 infinity-corrected (∞) objectives require a 14111 adapter.*  
*Microscopes with Leica 25 mm thread/45 mm parfocal length objectives require a 14115 adapter.*  
*Microscopes with 160 mm tube length objectives do not require an adapter.*



**Objectives, Jenco, Meiji, National, Nikon, Olympus**

160 mm tube length, DIN

Achromats, Motic/National

4X, N.A. 0.10. <b>704-160</b> .....	\$30.00
10X, N.A. 0.25. <b>710-160</b> .....	\$31.00
40X, N.A. 0.65. <b>740-160</b> .....	\$40.00
100X Oil, N.A. 1.25. <b>799-160</b> .....	\$65.00

Super High Contrast Achromats, National

4X, N.A. 0.10. <b>704-160ASC</b> .....	\$26.00
10X, N.A. 0.25. <b>710-160ASC</b> .....	\$34.00
40X, N.A. 0.65. <b>740-160ASC</b> .....	\$62.00
100X Oil, N.A. 1.25. <b>799-160ASC</b> .....	\$95.00

Semi-Plan Achromats, National (unless noted)

2X, N.A. 0.05, Accu-Scope. <b>3193SP-00</b> .....	\$70.00
4X, N.A. 0.10. <b>704-160SP</b> .....	\$25.00
10X, N.A. 0.25. <b>710-160SP</b> .....	\$48.00
40X, N.A. 0.65. <b>740-160SP</b> .....	\$95.00
50X Oil, N.A. 0.95, Meiji. <b>MA694</b> .....	\$280.00
60X, N.A. 0.85. <b>MA645</b> .....	\$300.00
100X Oil, N.A. 1.25. <b>799-160SP</b> .....	\$150.00

Plan Achromats, National (unless noted)

2.5X, N.A. 0.08, Meiji. <b>MA432</b> .....	\$270.00
4X, N.A. 0.10. <b>704-160P</b> .....	\$31.00
10X, N.A. 0.25. <b>710-160P</b> .....	\$52.00
20X, N.A. 0.40. <b>720-160P</b> .....	\$115.00
40X, N.A. 0.65. <b>740-160P</b> .....	\$99.00
50X Oil, N.A. 0.95, Accu-Scope. <b>3190PL-00</b> .....	\$315.00
60X, N.A. 0.85. <b>760-160P</b> .....	\$185.00
100X Oil, N.A. 1.25. <b>799-160P</b> .....	\$160.00