

# SLIDE & CASSETTE PROCESSING

## Microscope Slides

Slides, Economy, 75 x 25 x 1 mm, 72/bx, 20 bx/cs	<u>1 bx</u>	<u>20 bx</u>
Plain, ground edges. 9101 .....	\$3.40	65.00
Frosted (sandblasted), beveled. 8205 .....	\$5.05	95.00
Slides, premium, 75 x 25 x 1 mm, 72/bx, 20 bx/cs	<u>1 bx</u>	<u>20 bx</u>
Frosted (sandblasted), Erie. 2951 .....	\$10.00	168.00
Superfrost® (white), "SPECIMEN", Erie. 4951 .....	\$13.65	229.00

Superfrost® Plus

Concavity, 2-Well

Diagnostic Multi-Well



### Slides, specialty

#### Adhesion, 75 x 25 x 1 mm

Recommended for adhering tissue sections and cell preparations.

Starfrost® (silanated), 72/bx, 20 bx/cs,	<u>20 bx</u>
White end, Germany. 7255 .....	\$263.00
Superfrost® Plus (silanated), 72/bx, 20 bx/cs	<u>1 bx</u> <u>20 bx</u>
White end, Erie. 4951+ .....	\$29.15   488.00
<i>Preferred slide for Ventana instrumentation</i>	
Colorfrost® Plus (silanated), 72/bx, 20 bx/cs	<u>20 bx</u>
Blue end, Erie. 5951+ .....	\$499.00
Pink end, Erie. 7951+ .....	\$499.00

#### Cytogenetics, 75 x 25 x 1 mm

These slides are manufactured to provide a uniform chromosome spread.

Superfrost® CyGen, Erie. 4951-7140-006 .....	<u>20 bx</u>
	\$325.00

#### Concavity

75 x 25 x 1.2 mm, economy, 12/box	<u>ea</u>	<u>1 bx</u>
1-Well, 15 mm dia., 0.5 mm deep. 7103 .....	\$0.60	5.80
2-Well, 15 mm dia., 0.5 mm deep. 7104 .....	\$0.70	6.75
76 x 26 x 1.5 mm, Erie, 12/box	<u>1 bx</u>	
1-Well, 18 mm dia., 0.8 mm deep. 1518 .....	\$38.80	
2-Well, 18 mm dia., 0.8 mm deep. 1519 .....	\$50.60	
76 x 26 x 3 mm, Erie, 12/box	<u>1 bx</u>	
1-Well, 18 mm dia., 1.75 mm deep. 1526 .....	\$58.30	
2-Well, 18 mm dia., 1.75 mm deep. 1527 .....	\$70.85	

#### Diagnostic multi-well, 75 x 25 x 1 mm, 72/bx (10 bx min)

Cel-Line, white end, Teflon® printed .....	<u>1 bx</u>
	\$39.30
<i>(Please call regarding availability of configurations not listed. Specify color of mask, no. and diameter of wells.)</i>	
2-Well, 7 mm. 30-976H-BLUE	
8-Well, 5 mm. 30-962H-GREEN	
8-Well, 7 mm. 30-154H-BLUE	
10-Well, 6 mm. 30-225H-BLACK	
10-Well, 7 mm w/Numbers. 30-370H-BLUE	
18-Well, 5 mm. 30-943H-BLUE	

#### Fused Quartz, NSG Type OZ®, 3" x 1" x 1 mm

QMS3X1X1, each .....	\$20.00
----------------------	---------

#### Fused Silica, Corning 7980, 3" x 1" x 1 mm

FSMS3X1X1, each .....	\$29.00
-----------------------	---------

3 x 2", economy, 144/box. 6101 .....	\$15.10
--------------------------------------	---------

## Coverslips (Cover Glasses)

### Coverslips, premium, borosilicate glass, 1 oz/pkg

#### #0 Thickness (0.08-0.13 mm)

##### Squares & Rectangles, Erie

22 x 22 mm, ~240 pcs/oz, Cs/10 oz. 0-2222 .....	\$136.00
22 x 40 mm, ~132 pcs/oz, Cs/10 oz. 0-2240 .....	\$132.00

#### #0 Thickness (0.08-0.12 mm)

##### Circles, Germany\*

25 mm, ~222 pcs/oz, Cs/10 oz. 0-25 .....	\$243.00
--	----------

#### #1 Thickness (0.13-0.17 mm)

##### Squares & Rectangles, Erie

18 x 18 mm, ~242 pcs/oz. 1-1818 .....	—	55.70
22 x 22 mm, ~156 pcs/oz. 1-2222 .....	\$6.45	52.40
22 x 30 mm, ~116 pcs/oz. 1-2230 .....	—	52.40
22 x 50 mm, ~70 pcs/oz. 1-2250 .....	\$5.90	48.20
24 x 30 mm, ~108 pcs/oz. 1-2430 .....	—	52.40
24 x 50 mm, ~65 pcs/oz. 1-2450 .....	\$5.90	48.20
24 x 60 mm, ~54 pcs/oz. 1-2460 .....	—	48.20

#### #1 Thickness (0.13-0.17 mm)

##### Circles, Germany\*

12 mm, ~663 pcs/oz, Cs/5 oz. 1-12 .....	\$264.00
18 mm, ~296 pcs/oz, Cs/5 oz. 1-18 .....	\$127.00

#### #1.5 Thickness (0.16-0.19 mm). Recommended for most users

##### Squares & Rectangles, Erie

18 x 18 mm, ~209 pcs/oz. 1.5-1818 .....	—	55.70
22 x 22 mm, ~138 pcs/oz. 1.5-2222 .....	\$6.45	52.40
22 x 30 mm, ~102 pcs/oz. 1.5-2230 .....	—	52.40
22 x 40 mm, ~76 pcs/oz. 1.5-2240 .....	—	48.20
22 x 50 mm, ~61 pcs/oz. 1.5-2250 .....	\$5.90	48.20
22 x 60 mm, ~51 pcs/oz. 1.5-2260 .....	—	48.20
24 x 30 mm, ~94 pcs/oz. 1.5-2430 .....	—	52.40
24 x 40 mm, ~70 pcs/oz. 1.5-2440 .....	\$5.90	48.20
24 x 50 mm, ~56 pcs/oz. 1.5-2450 .....	\$5.90	48.20
24 x 60 mm, ~47 pcs/oz. 1.5-2460 .....	\$5.90	48.20

#### #1.5 Thickness (0.16-0.19 mm). Recommended for most users

##### Circles, Germany\*

12 mm, ~552 pcs/oz. 1.5-12 .....	\$238.00	—
22 mm, ~164 pcs/oz. 1.5-22 .....	—	\$238.00
25 mm, ~127 pcs/oz. 1.5-25 .....	—	\$257.00

##### Circles, Erie

40 mm, ~54 pcs/oz. 1.5-40 .....	\$26.20	211.00
---------------------------------	---------	--------

40 mm required for Biotech's FCS2 chamber.

\*German glass is preferred for cell growth/adhesion

Continued

## Coverslips, specialty

### Fused Quartz

0.15-0.25 mm Thick, NSG Type OZ®	
1" Square, 1 each. QCS1X1 .....	\$19.00
1" Circle, 1 each. QCS1X0 .....	\$19.00
0.16 ± 0.01 mm Thick, A1	
25 mm Square, 1 each. R525025 .....	\$32.40
25 mm Circle, 1 each. R525000 .....	\$23.40

### Fused Silica


0.15-0.25 mm Thick, Corning 7980	
1" Square, 1 each. FSCS1X1 .....	\$21.50
1" Circle, 1 each. FSCS1X0 .....	\$21.50
0.16 ± 0.01 mm Thick, S1-UV	
25 mm Square, 1 each. R425025 .....	\$32.00
25 mm Circle, 1 each. R425000 .....	\$28.80

### Superslips®, #1 thickness (0.13-0.17 mm), Erie

*Recommended for automatic coverslipppers and humid environments.*

	<u>10 oz</u>
24 x 40 mm .....	\$62.20
24 x 50 mm .....	\$62.20
24 x 60 mm .....	\$62.20


## Accessories



	<u>ea</u>	<u>cs/48</u>
Bibulous Paper, 10 x 15 cm, Bk/50 Sheets. PB46L .....	\$3.00	130.00



Spade Tip



Square

### Coverslip Forceps

Spade Tips, angled	
Stainless Steel, 11 cm. 1243 .....	\$2.25
Teflon® Coated, 11 cm, Bel-Art. 379420000 .....	\$14.40
Square Point Tips, straight	
Stainless Steel, 12 cm. 1244 .....	\$2.25
Teflon® Coated, 11 cm, Bel-Art. 379430000 .....	\$16.15



Mount-Quick

### Mounting Media

#### Aqueous-Based

Fluoromount-G™, 25 mL, SBA. 100-01 .....	\$70.00
Reduce fluorochrome quenching on slides stained with an aqueous final step.	
Mount-Quick, 30 mL in Dropper Bottle, Kiyota. K-750 .....	\$18.00
Suitable for fluorescence microscopy.	

#### Solvent-Based (Organic)

Mount-Quick, 30 cps, 30 mL in Tube, Kiyota. K-700 .....	\$12.00
Takes only 10 minutes to dry.	

Sealant, Hard as Nails®, 13.3 mL. 4400-66 .....	\$3.00
Seal edges of coverslip with clear nail polish.	



25-Slide  
Squeeze-Release Cover

## PLASTIC SLIDE BOXES



100-Slide  
Steel-Hinged Cover



50-Slide  
Hinged Cover

### Slide Boxes, hold 75 x 25 mm (3 x 1") slides

#### 25-Slide

Squeeze-Release cover, cork bottom, 9.5 x 8.5 x 3 cm	
Blue, ABS plastic, China. 97-0026 .....	\$1.75
Green, ABS plastic, China. 97-0027 .....	\$1.75
Red, ABS plastic, China. 97-0028 .....	\$1.75
White, ABS plastic, China. 97-0025 .....	\$1.75
Yellow, ABS plastic, China. 97-0029 .....	\$1.75
One-Piece construction, beige, USA. PP-25 .....	\$3.40

Resistant to freezer damage and light penetration.

UV protected polypropylene. Dimensions: 12 x 9.5 x 4 cm.

Friction-Fit cover, cork bottom, 14.1 x 8.8 x 3.5 cm, China

White, ABS plastic. HS15989E .....	\$5.30
------------------------------------	--------

#### 50-Slide, plastic, hinged cover, cork bottom, 20 x 9.5 x 3.5 cm

Blue, ABS plastic, China. 97-0056 .....	\$5.00
Red, ABS plastic, China. 97-0058 .....	\$5.00
White, ABS plastic, China. 97-0055 .....	\$5.00

#### 100-Slide

One-Piece construction, unlined	
Blue, polypropylene. HS15991A .....	\$5.00

Steel-Hinged cover, Cork Bottom	<u>ea</u>	<u>cs/48</u>
Gray, polystyrene, USA. PS-100 .....	\$9.00	388.00

PS-100 high-impact slide box offers superior resistance to freezer damage and light penetration. UV protected.

Blue, ABS plastic, China. HS15994A .....	\$6.80
Green, ABS plastic, China. HS15994B .....	\$6.80
Red, ABS plastic, China. HS15994C .....	\$6.80
White, ABS plastic, China. HS15994E .....	\$6.80
Yellow, ABS plastic, China. HS15994D .....	\$6.80

#### Steel-Hinged cover, Foam Bottom

White, ABS plastic, China. HS15994R .....	\$5.80
---	--------

Wood Frame Covered w/Black Paper .....	\$12.30
--	---------

A plastic, numbered insert holds slides, USA.

Metal Slide Cabinet



Wood Slide Cabinets

**Slide Cabinets**

Metal (steel), Paragon

(\*Specify Green, Gray, or Tan for cabinet and/or base)

6-Drawer cabinet, stackable. C101\_\* .....\$114.00

Holds 4500 75 x 25 slides in 6 divided drawers.

Memorandum Markers, paperboard, Pk/100. C108-M \$6.00

Compressor Follow Block, rubber, each. C109 .....\$1.15

12 Cabinet System (includes base). C112B\_\* .....\$1,372.00

Base, raises cabinets 4 3/4" off floor. C103\_\* .....\$40.00

Wood (oak) with Aluminum Trays, Eberbach

Holds 75 x 25 mm slides

500-Slide Capacity. 2800 .....\$523.00

1000-Slide Capacity. 2820 .....\$984.00

2000-Slide Capacity. 2830 .....\$1,823.00

3000-Slide Capacity. 2840 .....\$2,447.00

Slide Tray only, holds 20 slides. 2870 .....\$20.60

Fiberboard File



Index Card File



**Slide Files**


Fiberboard, hinged lid, Cargille. 73010 .....\$31.50

Holds 480 75 x 25 mm slides, or 240 75 x 51 mm slides in indexed compartments. Dimensions: 24.1 x 15.2 x 8.3 cm.

Index Cards with File Box. HS9901B .....\$12.30

12 Tabbed index cards hold 4 slides each in plastic pockets.


Index Cards only, Pk/12. HS9901A .....\$11.00



Holds up to 200 slides vertically in 100 individually numbered slots. Dimensions: 8.2 x 24.5 x 8.6 cm (H).

Polystyrene, tinted hinged cover, Cs/10. M700-100\_\* .....\$99.00

(\*Specify Blue, Green, Pink, or Yellow base color)



File Box

Slide Sleeves, 10 gauge vinyl

4 x 5 1/2" Card Holder Plus 4 Slide Pockets. SL-3-264 .....\$0.45

File Box, 4 x 6 x 11 1/2", fiberboard. 40589 .....\$16.80

5 x 8" Card Holder Plus 5 Slide Pockets. SL-3-108 .....\$0.55

File Box, 6 x 9 x 11 1/2", fiberboard. 40591 .....\$26.25

8 1/2 x 11" Page Holds 18 Slides. SL-3-182 .....\$0.75

Cardboard Slide Folder



Plastic Slide Folder

**Slide Folders, 20-slide**

Cardboard, USA

Black Tape. SF20-BLACK .....\$6.90 432.00

Canary Yellow Tape. SF20-YELLOW .....\$6.90 432.00

Fire Engine Red Tape. SF20-RED .....\$6.90 432.00

Forest Green Tape. SF20-GREEN .....\$6.90 432.00

Navy Blue Tape. SF20-BLUE .....\$6.90 432.00

*Folders made in the USA resist heat-induced warping* .....ea

Cardboard, Black Tape, China. HS9902 .....\$4.50

Plastic, gray tray w/2 clear doors, USA. 62680 .....\$9.15

Transparent covers hold visible slides securely in place, but can fold flat against back of tray when opened.

Slides are slightly elevated for easy retrieval.



Hemostat Slide Forceps



Kirkbridge Slide Forceps

**Slide Forceps**

Minimize contamination during slide handling.

Polypropylene. Will not scratch; corrosion and heat resistant.


Hemostat, 12.5 cm, Sweden. HEM-300.00 .....\$1.25

Straight, serrated handles and tips, 13 cm. 97-3901 .....\$0.60

Wire, Kirkbridge, nickel-plated, 12.5 cm long. 7-G36 .....\$1.95



12-Slide



52-Slide



91-Slide


**Slide Holders, perfect for slide drying or staining**

12-Slide (horizontally positioned), polystyrene. 445570000 .....\$5.95

52-Slide (approx.), HDPE. L500-C .....\$27.20

91-Slide (approx.), HDPE. L500-D .....\$36.35

Cardboard Mailer  
2-Slide



Plastic Mailer  
2-Slide



Plastic Mailer  
5-Slide



**Slide Mailers**

Cardboard, top corner notched	ea	cs/2000
1-Slide. MSL-01 .....	\$0.14	231.00
	ea	cs/1000
2-Slide. MSL-02 .....	\$0.23	189.00
Plastic	ea	cs/1000
1-Slide. 880052 .....	\$0.13	102.00
	ea	cs/500
2-Slide. 880041 .....	\$0.30	116.00
	ea	cs/200
4-Slide, screw-on cap. 950240 .....	—	90.00
	ea	cs/200
5-Slide, end opening, USA. 62706 .....	\$0.45	78.00

62706 mailers can also be used as staining jars.

**Slide Markers**



Pens,		
Pigment Liner, 0.5 mm stroke, Staedtler	ea	bx/10
Black. 308-05-9 .....	\$3.50	\$29.00

Write a 0.4 mm stroke on glass slides with this permanent marker.  
Ideal for cytologists and pathologists.



Pilot®, permanent, ultra fine point	ea	bx/12
Black. 44102 .....	\$2.22	21.60
Blue. 44103 .....	\$2.22	21.60
Green. 44105 .....	\$2.22	21.60
Red. 44104 .....	\$2.22	21.60

Ink resists xylene, alcohol, acetone and formalin.  
Remove with acetic acid, chloroform, or ethanol/xylene mixture.



Securline® Marker II®, permanent, PDC	ea	bx/10
Black, 0.6 mm tip. 1451 .....	\$2.50	20.30



Confine immunohistochemical reagents within a light blue, hydrophobic barrier, or separate multiple specimens on a slide/dish. Resists alcohol and acetone, but removes with xylene.

Super <sup>HT</sup> PAP (Peroxidase Anti-Peroxidase), hydrophobic		
Regular Tip (4 mm tapered nib). K-600 .....	\$54.00	
Fine Tip (2.5 mm tapered nib). K-650 .....	\$41.00	



Wet Surface. 133820000..... \$23.50



Pen Holder, Captor™, Bel-Art, Pk/3. 133730000..... \$15.35



Scribes, retractable		
Carbide-Tipped (resharpenable). RSC .....	\$15.30	
Replacement Carbide Tip. RCT-1 .....	\$6.20	
Diamond-Tipped. RSD .....	\$19.10	
Replacement Diamond Tip. RDT-1 .....	\$12.90	

Aluminum  
Slide Tray



**Slide Trays**

Aluminum, Eberbach. 2870 .....	\$20.60
Plastic, white, 20-slide, Bel-Art. 445810000 .....	\$9.45
Wood, 50-slide .....	\$42.00

Plastic Slide Tray



**Slide Warmers**

Small, 26 x 18 cm (10 x 7"). XH-2002 .....	\$225.00
Large, 63 x 20 cm (25 x 8"). XH-2001 .....	\$299.00


## Immersion Oil



Cargille Immersion Oil

<b>Type A, Low Viscosity (150 cSt). 16482. Short W.D. optics</b>	
7.5 mL in Dispensing Bottle. 18101 .....	\$2.90
20 mL in Dispensing Bottle. 18102 .....	\$7.75
120 mL. 18103 .....	\$13.80
480 mL. 18104 .....	\$30.25
<b>Type B, High Viscosity (1250 cSt). 16484. Long W.D. optics</b>	
7.5 mL in Dispensing Bottle. 18105 .....	\$2.90
20 mL in Dispensing Bottle. 18106 .....	\$7.75
120 mL. 18107 .....	\$13.80
480 mL. 18108 .....	\$30.25
<b>Type NVH, Very High Viscosity (21,000 cSt). 16485</b>	
120 mL. 18112 .....	\$18.45
480 mL. 18113 .....	\$42.40
<b>Type DF, Low Fluorescence (330 cSt). 16242</b>	
7.5 mL in Dispensing Bottle. 18115 .....	\$4.65
20 mL in Dispensing Bottle. 18116 .....	\$8.70
120 mL. 18117 .....	\$21.80
480 mL. 18118 .....	\$51.40
<b>Type FF, Fluorescence Free (170 cSt). 16212</b>	
7.5 mL in Dispensing Bottle. 18121 .....	\$4.45
20 mL in Dispensing Bottle. 18122 .....	\$8.60
120 mL. 18123 .....	\$20.85
480 mL. 18124 .....	\$49.00
<b>Type 37, 37°C (1250 cSt @ 37°C). 16237</b>	
20 mL in Dispensing Bottle. 18125 .....	\$8.30
120 mL. 18126 .....	\$18.45
480 mL. 18127 .....	\$42.40

A long nozzle, textured sides, and slim profile make this chemically inert, high-density polyethylene (HDPE) container ideal for dispensing immersion oil.



Dispensing Bottle, HDPE, 20 mL (Empty). 18145 .....\$5.80



Store immersion oil or mounting media in this borosilicate glass bottle with an applicator rod and ground glass cap.


Balsam Bottle, glass, 45 mL (Empty), USA. 616 .....\$14.80

## Staining

Glass staining jar with ground glass lid holds 8 coverslips, 18-22 mm, back-to-back. Jar is 45 mm wide and tall.




Coverslip Staining Jar. 200-011 .....\$11.35



<b>Staining/Moisture Chambers, 10-slide, USA</b>	<u>ea</u>	<u>pk/3</u>
Clear, polystyrene. 240-9000-010 .....	\$47.50	133.00
Black, polystyrene. 240-9020-Z10 .....	\$48.80	137.00



<b>Staining Dishes, hold slides horizontally</b>	
Glass (flint), 20-slide (up to 76 x 26 mm), back-to-back, USA	
Set (grooved dish & lid). 672 .....	\$12.35
Bottom only. 672-BTM .....	\$8.65
Lid only. 672-LID .....	\$4.80



<b>Plastic</b>	
5-Slide, ABS plastic. HS15989 .....	\$3.70
10-Slide, ABS plastic. HS15987 .....	\$4.00



Plastic Set

Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).

25-Slide, POM (polyoxymethylene) plastic	
Set (dish, rack & lid). 0990-99 .....	\$48.75



674



442081000

675



<b>Staining Jars, hold slides (up to 76 x 26 mm) vertically</b>	
Coplin, holds 10 slides back-to-back	
Glass w/Ground Glass Cover, USA. 673 .....	\$10.90
Replacement Glass Cover. 673-LID .....	\$3.95
Glass w/Screw-On Phenolic Cap, USA. 674 .....	\$10.90
Replacement Phenolic Cap, black. 674-LID .....	\$4.10
Polypropylene w/Screw-On Cap, Bel-Art. 442081000 .....	\$7.10
Plastic jar is unbreakable, microwave proof, stain resistant, and chemically inert.	
<b>Wide Top</b>	
Glass, holds 16 slides back-to-back	
Set (grooved jar & lid), USA. 675 .....	\$14.65
Lid only. 675-LID .....	\$4.30

*Continued*



12-Slide  
Move slides without forceps between jars linked together in a workstation. Swing out the rack's integral lid for easy accessibility to slides. Stable, wide-based jar requires only 80 mL of reagent. Polyoxymethylene (POM) resists alcohol and xylene; can be microwaved or frozen.

Plastic (POM)

12-Slide

Jar only

	ea	cs/6
Blue. M900-12B	\$11.35	63.15
Green. M900-12G	\$11.35	63.15
Pink. M900-12P	\$11.35	63.15
White. M900-12W	\$11.35	63.15
Yellow. M900-12Y	\$11.35	63.15

Rack/Lid only, gray. M905-12DGY .....\$12.60 70.20

25-Slide, light green

Jar only. 1030 .....\$10.80

Rack only

Brass Handle. 1040 .....\$15.20

Plastic Handle (microwaveable). 1042.....\$17.30

Tissue (Embedding) & Biopsy Cassettes




Base Molds, disposable, Simport

15 x 15 x 5, Cs/1000. M475-2	\$134.00
24 x 24 x 5, Cs/1000. M475-3	\$134.00
30 x 24 x 5, Cs/1000. M475-4	\$134.00
37 x 24 x 5, Cs/1000. M475-5	\$134.00




Biopsy Bags (Tea Bags), paper, 2 1/4 x 2 1/4", Pk/500. 56100.....\$39.00

Made of 60 PPI cellular polyester urethane foam.



Biopsy Pads

Black Foam, 34.3 x 25.4 mm, Bg/1000. 56150	\$37.50
Blue Foam, 30.2 x 25.4 x 2 mm, Cs/10,000. M476-1	\$191.00



Molded from a high-density, solvent resistant polymer. Textured writing surface on front and sides.

Cassettes, Detached Lids (packed separately)

Biopsy, small square pores, economy, Pk/500.....\$39.50

Blue. BC-0109B White. BC-0109W

Green. BC-0109G Yellow. BC-0109Y

Tissue, slotted pores

Economy, Pk/500 .....\$39.50

Blue. EC-0109B Pink. EC-0109P

Green. EC-0109G White. EC-0109W

Lilac. EC-0109L Yellow. EC-0109Y

Uni-Cassettes™, premium, Pk/500 .....\$57.50

Aqua. 6059 Lime. 6062

Blue. 6054 Peach. 6064

Burnt Orange. 6056 Pink. 6052

Chartreuse. 6061 Red. 6060

Gold. 6058 Rose. 6063


Green. 6053 Tan. 6055

Lavender. 6065 White. 6050

Lilac. 6057 Yellow. 6051



Histosette I



Micromesh

Cassettes, Attached Lids, acetyl polymer, premium

Biopsy

Histosette™ I, 1 mm square pores, Pk/1500 .....\$141.00 381.00

Aqua. M491-12 Peach. M491-7

Blue. M491-6 Pink. M491-3

Green. M491-4 Tan. M491-8

Lilac. M491-10 White. M491-2

Orange. M491-11 Yellow. M491-5

Micromesh™, 0.38mm pores, Pk/1000 .....\$249.00 672.00

Aqua. M507-12 Peach. M507-7

Blue. M507-6 Pink. M507-3

Green. M507-4 Tan. M507-8


Lilac. M507-10 White. M507-2

Orange. M507-11 Yellow. M507-5

Continued



Histosette I



Macrosette


Tissue, slotted pores  
 Histosette™ I, snap-apart lids, Pk/1500 ..... 1 pk 3 pk  
 \$115.00 312.00  
 Aqua. M490-12 Peach. M490-7  
 Blue. M490-6 Pink. M490-3  
 Green. M490-4 Tan. M490-8  
 Lilac. M490-10 White. M490-2  
 Orange. M490-11 Yellow. M490-5

Uniset™, hinged lids, Pk/1500..... 1 pk 3 pk  
 \$124.00 336.00

(Alternative to Sakura's Uni-Cassette®)

Aqua. M505-12 Peach. M505-7  
 Blue. M505-6 Pink. M505-3  
 Green. M505-4 Tan. M505-8  
 Lilac. M505-10 White. M505-2  
 Orange. M505-11 Yellow. M505-5

Macrosette™, 2X normal depth, white, Pk/750. M512.....\$173.00



**Cassettes, Quickload™ stacks, acetyl polymer**  
 (Use in Leica and Sakura labelers)

Biopsy  
 Uniset™, hinged lids, Pk/1000..... 1 pk 3 pk  
 \$107.00 290.00  
 Aqua. M506-12T Peach. M506-7T  
 Blue. M506-6T Pink. M506-3T  
 Green. M506-4T Tan. M506-8T  
 Lilac. M506-10T White. M506-2T  
 Orange. M506-11T Yellow. M506-5T



M480-4

**Cassettes, without Lids** (compatible with all available metal lids)  
 Tissue, acetyl polymer, premium  
 Capsettes™, Pk/500.....\$39.60  
 Aqua. 6039 Lime. 6042  
 Blue. 6034 Peach. 6044  
 Burnt Orange. 6036 Pink. 6032  
 Chartreuse. 6041 Red. 6040  
 Gold. 6038 Rose. 6043  
 Green. 6033 Tan. 6035  
 Lavender. 6045 White. 6030  
 Lilac. 6037 Yellow. 6031  
 Simport, Pk/1500  
 Green, M480-4.....\$89.50



Store 1500 cassettes or 990 embedding rings in 6 stackable drawers.  
 Drawer dimensions: 40.5 x 23 x 5.1 cm.


**Cassette Cabinet, polystyrene, 6-drawer. M495-6** ..... \$147.00



**Freeze Spray, 10 oz Can. 403A-285G** ..... \$11.00

**Paraffin Scrapers** see Scrapers in LAB, SURGICAL, & DENTAL INSTRUMENTS section

Store 100-1000 blocks in 100 indexed compartments.  
 Dimensions: 38.7 x 26.7 x 6.0 cm.



**Tissue Block File, fiberboard, Cargille. 72010**..... \$32.00

**Tissue Block Mailing Sleeve, Block Tunnel™** ..... \$2.25

**Equipment**



**Paraffin Dispenser, 2.5 gallon. XH-4002**..... \$758.00



**Tissue Floatation Bath. XH-1001** ..... \$182.00